

Product datasheet for **RN203214**

Rab1b (NM_001109979) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rab1b (NM_001109979) Rat Untagged Clone
Tag: Tag Free
Symbol: Rab1b
Synonyms: Rab1bl
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN203214 representing NM_001109979
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAACCCGAATATGACTACCTGTTTAAAGCTGCTTTTGATTGGTGACTCGGGCGTGGCAAGTCATGCC
TGCTTCTGCGTTTGCTGATGACACGTACACAGAGAGCTACATCAGCACCATTGGGGTGGACTTCAAGAT
TCGAACCATTGAAGTGGATGGCAAACCATCAAACACAGATTTGGGACACAGCTGGTCAGGAACGGTTC
AGGACCATCACTCCAGCTACTATCGGGTGCTCATGGCATCATTGTGGTGTATGACGCTACTGACCAGG
AATCCTACGCTAATGTGAAACAGTGGCTGCAGGAAATAGATCGCTACGCCAGTGAGAATGTCAATAAACT
GCTGGTAGGCAACAAGAGTGACCTCACCAAGAAGGTCGTGGACAATACCACAGCCAAGGAATTTGCA
GACTCTCTGGGTGTCCTTCTGGAGACAAGTGCCAAGAATGCCACCAATGTTGAACAGGCATTATGA
CAATGGCTGCAGAGATCAAAAAACGGATGGGGCCAGGAGCAGCATCTGGGGTGAACGGCCCAACCTGAA
GATCGACAGCACTCCTGTGAAATCTGCTAGTGGTGGCTGCTG**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
ACCN: NM_001109979
Insert Size: 606 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001109979.2, NP_001103449.1</u>
RefSeq Size:	1873 bp
RefSeq ORF:	606 bp
Locus ID:	100126191
Cytogenetics:	1q43