

Product datasheet for RN203286

Slc25a26 (NM_001108638) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Slc25a26 (NM_001108638) Rat Untagged Clone
Tag: Tag Free
Symbol: Slc25a26
Synonyms: SZGENESYMBOL2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN203286 representing NM_001108638
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGACGCGCCGGGCTTCACGGCCTCACTGGTGGCTGGAGGAGTAGCAGGTGCCGCTGTGGACTTGATAT
 TATTTCCTTGGATACTATTAAGACAAGGCTGCAGAGTCCCCAGGGATTTAACAGGCTGGCGGCTTTTCG
 TGGAAATATGCTGGTCCCTTCTGCAGCTATTGGATCCTTTCTAATGCTGCTGCGTTTTTTCTCACC
 TATGAATATGTAAATCTTTGCTACACACCGACTCAACTTCACACTTCAGACCCATGAAGCACATGCTGG
 CCGCCTCTACAGGAGAAGTGGTGCATGGGTGTTTCTGTAAGTGTGGCAAGTATGGACATATCTGTGGAG
 ACCACAGTTCACCTTAGGTCT**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
ACCN: NM_001108638
Insert Size: 375 bp
OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108638.1</u> , <u>NP_001102108.1</u>
RefSeq Size:	2228 bp
RefSeq ORF:	375 bp
Locus ID:	362403
Cytogenetics:	4q34