

## Product datasheet for **RN203307**

### **Actr1a (NM\_001106364) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Actr1a (NM\_001106364) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Actr1a  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN203307 representing NM\_001106364  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAGTCTACGATGTGATCGCCAACCAGCCTGTTGTGATCGACAATGGATCCGGTGTGATTAAGCTG  
GTTTTGCTGGTGATCAGATCCCCAAATACTGCTTTCCAACTATGTGGCAGACCCAAACATGTTTCGAGT  
CATGGCAGGAGCCCTCGAAGGTGACATCTTCATTGGCCCAAAGCTGAGGAGCACCGGGGCTGCTGTCA  
ATCCGCTACCCCATGGAACATGGCATCGTCAAGGACTGGAACGACATGGAGCGCATCTGGCAATATGTCT  
ATTCTAAGGACCAGCTGCAGACTTTCTCAGAGGAGCATCCTGTGCTCCTGACTGAGGCACCTTTAAACCC  
TCGGAAAAACCGGAGCGAGCGGCAGAAGTTTTCTCGAGACCTCAATGTGCCAGCTCTCTTCATCTCC  
ATGCAAGCTGTGCTTAGCCTTTATGCCACAGGCAGAACCACAGGTGTGGTGTGGATTCTGGGGATGGCG  
TCACCCATGCAGTCCCGATTTACGAAGGCTTCGCCATGCCTCACTCCATCATGCGCATTGACATTGCCGG  
CCGCGATGTCTCAGCTTCTCCGCTCTACCTACGAAAGGAGGGCTACGATTTCCACTCCTCCTCCGAG  
TTTGAGATTGTCAAGGCCATAAAGGAAAGAGCTTGCTACCTGTCCATAAACCCCAAGGATGAGACAC  
TGGAGACGGAGAAGGCTCAGTACTACCTGCCCGATGGCAGCACCATTGAGATTGGCCCTCCAGGTTCCG  
GGCCCTGAGCTGCTGTTCAAGCCTGACTTGATTGGCAGGAGAGTGAGGGTATCCATGAGGTGCTGGTG  
TTCGCTATCCAGAAGTCAGACATGGATTTACGGCGGACGCTCTTTTCCAACATTGTCTCTCAGGAGGTT  
CTACCCTCTTCAAAGTTTTGGGACAGGCTACTGAGTGAAGTGAAGAACTGGCTCCAAAAGATGTGAA  
GATCAGGATATCTGCACCCAGGAGAGACTGTATTCCACATGGATTGGAGGCTCTATCCTTGCCCTCCCTG  
GATACCTTTAAGAAGATGTGGTCTCCAAGAAGGAATATGAAGAAGATGGTGCCCGATCCATCCACAGGA  
AAACCTTCTAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001106364



[View online >](#)

<b>Insert Size:</b>	1131 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001106364.1</a></u> , <u><a href="#">NP_001099834.1</a></u>
<b>RefSeq Size:</b>	2767 bp
<b>RefSeq ORF:</b>	1131 bp
<b>Locus ID:</b>	294010
<b>UniProt ID:</b>	<u><a href="#">P85515</a></u>
<b>Cytogenetics:</b>	1q54
<b>Gene Summary:</b>	Component of a multi-subunit complex involved in microtubule based vesicle motility. It is associated with the centrosome (By similarity).[UniProtKB/Swiss-Prot Function]