

Product datasheet for RN203328

Olr44 (NM_001000125) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Olr44 (NM_001000125) Rat Untagged Clone
Tag: Tag Free
Symbol: Olr44
Synonyms: ratchr1-160300212-160299265_ORF
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN203328 representing NM_001000125
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCAACACACGAACCACAGCCATCAGAACCCAAGCTCCTTCCTGCTCATGGGGATTCCTGGCCTGGAGG
 CAGCCCACTTTGGATTGCATTTCTTTCTGTTCCATGTATGCCCTAGCAGTTCTTGCAACATGGCAGT
 CCTGCTGGTGGTGCCTCTGAGCCTTCACTGCACCAGCCCATGTACCTGTTCTGTGCATGCTGTCCACT
 ATTGACCTGATACTCTGCACCTCCACTGTCCCAAACTCCTCGCACTCTTTGGGCAAATGCAGCTGAGA
 TTGCCTTTGGGGCCTGTGCTGCCAGATGTTCTTTATCCATGGCTTCTCAGCTGTAGAATCTGGCATCCT
 GTTATCAATGGCCTTTGATCGATACTTGCCCATCTGCCGGCCACTGCACTATGGGTCACTGCTGTCCGCA
 GAGTCAGTGAGCAAGTTAGGGGCTGCTGCTCTGCTCCGTGGCCTGGGGCTCATGACCCCACTCACTTGCC
 TACTGGCAAGGCTGAGCTACTGTGGCCGTGTGGTAGCACACTCCTACTGTGAGCACATGGCTGTGGTGAA
 GCTGGCTGTGGGGGACCCAACCAACAATATATGGCATCACAGCAGCCACACTAGTTGTAGGAACT
 GACTCCATTTGTATTGCTGTATCCTATGCACTCATCTCCGGGCTGTGTTGGCCTCTCCTCAAAGAAG
 CTAGGGCAAAGACCTTTGGCACTTGTGGCTCCACCTGGGCGTCATCCTTCTTCTACACACCAGGGCT
 CTTCTCATTCTATACGCAGCGTTTCGGTCAACACGTACCCAGGCACGTCCACATCCTCCTGGCTGACCTC
 TACCTGGTTGTACCACCCATGCTCAACCCCATCATCTATGGGATGAAAACCAAGCAGATCAGAGACGGAG
 CACTCCGACTACTCAAGAGGGGCCCTGCACAATCA**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001000125



[View online »](#)

Insert Size:	948 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001000125.1, NP_001000125.1</u>
RefSeq Size:	948 bp
RefSeq ORF:	948 bp
Locus ID:	293193
Cytogenetics:	1q32
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]