

Product datasheet for RN203330

Olr820 (NM_001000974) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olr820 (NM_001000974) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Olr820
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN203330 representing NM_001000974 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAAGAAAATCAGACCATGGTCACAGAGTTCATCCTGCTGGGATTCTGTCTTGCCCTAGGATTACAC
TAGTTCTTTTCTACTTTTCTCTCTCTTCTATACGTTTACCATACTGGGCAATGGGACTATCCTTGCACT
GATCTGCCTAGACTCCAGACTCCACACTCCCATGTACTTCTTCTGTCCACCTGGCCATTGTCGATATG
GCCTATGCCTGCAACACAGTGCCACAAACACTGGTAAACCTCTTGGATGAGACCAAGTCCATCTCCTTCA
CTGGATGCATAACACAGACCTTTCTCTTCTTGGCTTTTGCACACAGAATGTGTCTCCTTGTCTGAT
GTCCTATGACCGGTATGTAGCTATCTGCCATCCACTACACTACACTGTCATTATGAATTGGAGAGTGTGT
ACCATCTTGTTGCTATTTCTGGATATTTAGCTTCTCTTGGCTCTGGTCCATATAATCCTCATCCTGA
GGCTTCCCTTCTGTGGACCTCATGAAATCAATCACTTCTTCTGTGAAATCCTGTCTGTCTCAAGCTGGC
CTGTGCTGACACAACACTCAATCAAGTTGTCATCTTGCAGCCTGTGTATTGATCTTAGTGGGGCCCTG
TGCTTTGTTCTAGTCTCCTATACATGCATCCTGATGGCCATCCTAAAGATCCAGTCAGGGGAGGGCCGCA
GAAAGGCCTTCTCCACCTGTTCTCCACCTCTGTGTGGTAGGGCTCTTCTTGGCAGTGCCATTGTCAT
GTACATGGCCCCAAGTCAGAGCACCTGAGGAACAGCAGAAGATCCTTTTCTTGTTTTACAGTTTTTTC
AACCCCATGCTGAACCCCTGATCTACAGCCTGAGGAATGCTGAGGTCAAAGGTGCCCTCTGGAGGTAC
TGTGCAAAGAAAGTCATTCTTGTGGTGTGGTGTTCACACCACAAGGCT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001000974
Insert Size:	963 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001000974.1, NP_001000974.1</u>
RefSeq Size:	963 bp
RefSeq ORF:	963 bp
Locus ID:	405329
Cytogenetics:	4q24
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]