

## Product datasheet for **RN203341**

### **Rab27b (NM\_053459) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Rab27b (NM\_053459) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Rab27b  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN203341 representing NM\_053459  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGACTGATGGAGACTATGACTACCTGATTAACCTCTGGCCCTTGGAGACTCAGGAGTAGGAAGACAA  
CATTTCTTACCAGATACACAGACAATAAATTCAATCCAAATTCACACTACAGTAGGAATTGATTTTCG  
GGAAAAACGTGTGGTTTATGACACACAAGGAGCAGATGGATCTTCAGGAAAGCCTTAAAGGTACATCTG  
CAGCTTTGGGACACTGCTGGACAAGAGCGGTTCCGGAGTCTCACCACCGCCTTCTCAGAGATGCCATGG  
GCTTCTTACTAATGTTTGACCTCACCAGTCAACAGAGTTTCTTGAATGTCAGAACTGGATGAGTCAACT  
GCAGGCAAATGCTTACTGTGAAAACCCAGATATAGTATTAATTGGCAACAAGGCAGACCTGCTGGACCAA  
AGGGAAGTCAACGAGCGACAAGCAGAGAGCTGGCTGAAAAATATGGCATAACCATACTTTGAAACAAGT  
CGGCGACAGGACAGAATGTAGAGAAGTCAGTGGAAACCTTCTGGACTTAATCATGAAGAGAATGGAGAA  
GTGTGTAGAGAAGACGAGGTTCCGGACACTGTCAATGGTGTGAATTCTGGGAAGGTCGATGGGGAGAAG  
CCAGCAGAGAAGAAATGTGCCTGCTAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_053459  
**Insert Size:** 657 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_053459.1</a></u> , <u><a href="#">NP_445911.1</a></u>
<b>RefSeq Size:</b>	890 bp
<b>RefSeq ORF:</b>	657 bp
<b>Locus ID:</b>	84590
<b>UniProt ID:</b>	<u><a href="#">Q99P74</a></u>
<b>Cytogenetics:</b>	18q12.2
<b>Gene Summary:</b>	member of the RAS oncogene family; may have small GTP binding activity [RGD, Feb 2006]