

Product datasheet for RN203365

Pus10 (NM_001025278) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pus10 (NM_001025278) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Pus10
Synonyms:	RGD1306402
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN203365 representing NM_001025278 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGTTCCCACTGAAGGAGGAGAACAGACGTGTGGCCAGCTGCTGCTCAGGTCTGGCACCTGCCAAGAT
 GCATCTTCAGATTCTGCGGCGTGGATCTTCCGGCCCTTACAAACGCCCATACAAGGAGTTGCTCGATGA
 TCTACAAACATTTCTGGAATCTGAGAACCGTGAGTTAATTCTGGAAGTCCAGACCCACCGTTGAAGAAG
 ATTCGTCTACACGAAGATGGGGTTGATACTTTGAGTGCAAATGGAGAGGAGGGCGTCTCTGTTACTGAAG
 ATGAAACATGGCGGCCACGAAGCCTTCAAATTGGAGAGTGTGCAATGTTTGCTTAGGGATTCTTCAGGG
 ATTCTGCGAGAAGGACTTTATCAGAGAGGTGTGCCAAAAGGTTGAAGCCTCTGGGTTTGAATTCACCAGC
 GTGGTTTTATCTGTGTCCTTCCACCGCAGCTGTCTGTGAGAGAGCATGCTGCATGGTTGCTGGTAAAC
 AGGAAATGGGAAAGCAGAGCCTGTCCCTGGGAAGAAATGATGTCGTTAGCTGAAGGAAGCCTACAAGTG
 GATAACTACCCCTGTTTTAGAGGAGCTGGGAGTCCCACTGATGGGAAGAGCTTGTTCAAGTGAGT
 GTGGTCTTTGCTACCCAGAAACACTTCAAGACTGCCATTTCTAGGTGAAGTTTGCCAGACTGTTTCA
 AACCAGCCAAAACAAACAGTCAGTATTCACCAGAATGGCCGTTTTGAAAGCGTTGAGTAGGATAAAAGA
 AGAAGACTTCCTCAGTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001025278
Insert Size:	789 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001025278.1, NP_001020449.1</u>
RefSeq Size:	1852 bp
RefSeq ORF:	789 bp
Locus ID:	305583
Cytogenetics:	14q22