

Product datasheet for **RN203400**

Rragd (NM_001106641) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rragd (NM_001106641) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Rragd
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN203400 representing NM_001106641 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGCCAGGTGCTGGGCAAGCCGACGCCGAGGGTGAGGACGGCGGCGACGACGAGGAGGAGGATGAGC
TGGTGGGGCTGGCAGGCTACGAGGATGGACCGGAGTCTCCGACGCCGAGCTGGACAGCGGACCGGAGGA
GGGAGTTCTGGACTTCAGCGACCCCTTCAGCACTGAGGTGAAGCCAGAATTCTGCTCATGGGTCTCAGG
AGGAGTGGCAAGTCGTCCATACAGAAAGTTGTCTTTCACAAAATGTCCCCGAGTGAGACGCTGTTCTCTGG
AGAGTACCAACAGGATCTGCCGGGAAGACGTCTCCAACAGCAGCTTTGTCAACTCCAGATTGGGACTT
CCCAGGGCAGATTGACTTTTTTGGCCCTACATTTGACTATGAGATGATCTCCGAGGGACAGGAGCACTG
ATCTTTGTCTAGACTCCAGGATGACTATATGGAGGCGCTGGCCAGGCTCCATCTCACTGTGACCAGGG
CCTACAAGGTGAACACTGATATCAACTTCGAGGTGTTTATTCATAAAGTGACGGTCTGTCGGATGATCA
CAAAATCGAAACCCAAAGAGACATTCACCAGAGGGCAAACGATGACCTTGCCGATGCTGGGTAGAGAAA
ATTCACCTGAGCTTTTATCTGACAAGCATATACGATCATTCAATATTCGAAGCCTTAGCAAAGTGGTTC
AGAACTGATCCCGCAGTCCCGACGCTGGAGAATCTGCTGAACATCTTTATCTCAAATCTGGAATTGA
AAAAGCATTTCTATTTGACGTGGTCAGTAAGATTTATATTGCAACCGACAGCACGCCGGTGGATATGCAG
ACGTATGAGCTGTGCTGCGACATGATTGATGTGGTGATTGACATCTCTTGATTTATGGCCTCAAAGAAG
ATGGAGCGGGCGCCCGTATGACAAGGACTCGACAGCCATTATAAACTGAATAATACAACGGTTCTTTA
TCTAAAAGAGGTACGAAGTTCTTGGCGCTCGTTTGCTTTGTGAGAGAGGAAAGCTTTGAAAGAAAAGGC
CTTATTGATTATAATTTCCACTGCTTCCGGAAGGCCATCCATGAAGTCTTCGAGGTGAGGATGAAGATGG
TGAAATCTCGGAAGGCCAGAGTCGGTTACCAAGAAGACTGGAGCTACTCCCAATGGGACTCCTCGGGT
GCTGCTG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001106641
Insert Size:	1200 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001106641.2, NP_001100111.1</u>
RefSeq Size:	1473 bp
RefSeq ORF:	1200 bp
Locus ID:	297960
Cytogenetics:	5q21