

Product datasheet for RN203433

Prdx5 (NM_053610) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Prdx5 (NM_053610) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Prdx5
Synonyms:	Aoeb166
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN203433 representing NM_053610 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGTCCAGCTGAGGTTTTGCGTCCTAGGCAGCATAGCCGGATCGGTGCTCCGTGCATCGGCTACTTGGACGTGCGTGGCAGGCAGAGCAGGCCGAAAGGAGCAGGTTGGGAGTGTGGTGGGGCCCGCAGCTTCAGCAGTGCCGCGGTGACTATGGCCCGATCAAGGTGGGAGACACCATTCCCTCAGTGGAGGTATTTGAAGGGGAACTGGAAGAAGGTGAACCTGGCAGAGCTGTCAAGGACAAGAAAGGTGTTTTGTTGGAGTCCCTGGGCATTTACACCTGGCTGTTCCAAGACCCATCTGCCTGGGTTTGTGGAGCAAGCCGGAGCTCTGAAGGCCAAGGGAGCACAAGTGGTGGCCTGTCTGAGTGTAAATGATGCCTTCGTGACTGCAGAGTGGGGTGGAGCCACAGGCAGAAAGCAAGGTTGAGCTCTGGCTGACCCACTGGAGCTTTGGAAAGGAGACAGATTTACTACTAGATGATTCTTTGGTGTCTCTCTTTGGGAATCGTCGGCTAAAAAGGTTCTCCATGGTGATAGACAAGGGCGTAGTAAAGGCACTGAACGTGGAGCCGGATGGCACAGGCCTCACCTGCAGCCTGGCCCCAACATCCTCTCACAACCT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_053610
Insert Size:	642 bp



[View online »](#)

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<p>NM_053610.1, NP_446062.1</p>
RefSeq Size:	<p>837 bp</p>
RefSeq ORF:	<p>642 bp</p>
Locus ID:	<p>113898</p>
UniProt ID:	<p>Q9R063</p>
Cytogenetics:	<p>1q43</p>
Gene Summary:	<p>member of the AhpC/TSA peroxiredoxin family with mitochondrial and peroxisomal sorting signals; may play an antioxidant protective role in various tissues [RGD, Feb 2006]</p>