

Product datasheet for RN203506

Gnai3 (NM_013106) Rat Untagged Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | Gnai3 (NM_013106) Rat Untagged Clone |
| Tag: | Tag Free |
| Symbol: | Gnai3 |
| Synonyms: | GIP3A |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >RN203506 representing NM_013106 Red=Cloning site Blue=ORF Orange=Stop codon |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGCTGCACGTTGAGCGCCGAGGACAAGGCGCGGTGGAGCGGAGTAAGATGATCGACCGCAACTTGC
GGGAGGACGGAGAGAAAGCGCCAAAGAAGTGAAGCTGCTGCTCGGCGCTGGAGAATCTGGTAAAAG
TACTATTGTGAAACAGATGAAAATCATTATGAGGATGGCTATCCGAGGACGAGTGAAGCAGTATAAA
GTTGTCGCTACAGCAATACCATTCACTCCATCATTGCAATCATAAGAGCCATGGGACGGTTGAAGATTG
ATTTTGGGAAGCTGCCAGAGCGGATGATGCCCGACAGTTATTTGTTTTAGCTGGCAGTGTGAGGAAGG
AGTCATGACTTCAGAAGTACGAGCGTGATTAACGTTTATGGCGAGATGGCGGGTGCAGGCATGCTTC
AGCAGGTCACAGGAATATCAGCTCAATGATTCTGCTTCATATTACCTAAATGATTTGGATAGAATATCCC
AGACCAACTACATTCCAAGTCAAGATGTTCTTCGGACGAGAGTGAAGACTACAGGCATTGTGGAGAC
CCACTTCACTTCAAGGAGCTCTACTTCAAAATGTTTGTAGTGTAGGTGGCCAAAGATCCGAACGAAAAAG
TGGATCCACTGTTTTGAGGGAGTGACAGCAATTATCTTTGTGTGGCTCTCAGTGATTACGACCTTGTTCT
TGGCTGAGGACGAGGAAATGAACCGAATGCATGAAAGCATGAAATGTTTGCAGCATTGTGAACAACAA
ATGGTTTACAGACACTTCAATCATTCTCTTCTTAATAAGAAAGACCTTTTTGAGGAAAAATAAAGAGG
AGTCCATTAACAATCTGTTATCCAGAATACACAGGTTCCAATACGTACGAAGAGGCAGCTGCTTACATCC
AGTGCCAGTTTGAAGATCTGAACCGAAGAAAGGACACCAAGGAGGTCTACACTCACTTTACCTGTGCCAC
TGACACCAAAAACGTGCAGTTTGTGTTTGTGCTGTTACAGATGTCATCATTAATAAACAACCTTAAAGGAA
TGTGGGCTTTACTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

| | |
|-------------------------------|--|
| ACCN: | NM_013106 |
| Insert Size: | 1065 bp |
| OTI Disclaimer: | <p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).</p> |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_013106.1</u> , <u>NP_037238.1</u> |
| RefSeq Size: | 3072 bp |
| RefSeq ORF: | 1065 bp |
| Locus ID: | 25643 |
| UniProt ID: | <u>P08753</u> |
| Cytogenetics: | 2q34 |
| Gene Summary: | has GTPase activity; involved in zymogen granule-plasma membrane fusion in pancreatic cells [RGD, Feb 2006] |