

Product datasheet for **RN203763**

Csrp1 (NM_017148) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Csrp1 (NM_017148) Rat Untagged Clone
Tag: Tag Free
Symbol: Csrp1
Synonyms: Csrp
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN203763 representing NM_017148
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCCAACTGGGGAGGAGGCAAGAAATGCGGGGTGCCAGAAGACAGTCTACTTTGCTGAGGAGGTCC
AGTGCGAGGGCAACAGCTTCCACAAATCCTGCTCCTGTGCATGGTTTGAAGAAGAATCTGGACAGCAC
CACTGTGGCAGTGCATGGAGAAGAGATCTATTGCAAGTCATGTTATGGCAAGAAGTATGGACCAAAAGGC
TACGGCTACGGGCAGGGTGCAGGCACGCTGAGCATGGACAAGGGGGAGTCTCTGGGCATCAAGCATGAGG
AAGCCCCGGACACAGGCTACCAACCAACCCCAATGCATCCAAGTTTGGCTCAGAAGATCGCGGGCTCTGA
GCGCTGTCCGCGATGCAGCCAGGCAGTCTATGCAGCGGAGAAGGTGATCGGTGCCGGGAAGTCTGGCAT
AAGTCCTGCTCCGATGTGCCAAATGTGGCAAAGGCCTCGAGTCGACCACCCTGGCAGACAAGGATGGTG
AGATCTATTGCAAAGGATGCTATGCCAAAACTTTGGCCCAAAGGTTTTGGCTTTGGACAAGGAGCTGG
CGCCTTGGTCCATTCAGAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_017148
Insert Size: 582 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_017148.2</u> , <u>NP_058844.1</u>
RefSeq Size:	1720 bp
RefSeq ORF:	582 bp
Locus ID:	29276
UniProt ID:	<u>P47875</u>
Cytogenetics:	13q13
Gene Summary:	cysteine rich protein that has homology to the delta opioid receptor; may play a role in neuronal development [RGD, Feb 2006]