

Product datasheet for RN203768

Upk3bl1 (NM_001109020) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Upk3bl1 (NM_001109020) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Upk3bl1
Synonyms:	Upk3bl
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN203768 representing NM_001109020 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGGGGCCCCACAGGAACAGTCTATGCTGCTGCTGCCCTGTTGGTGCTGGTGACTTGTGTCCAGCCCG
 GAACAGGTCTGGAGAGCATTAACTATGTCCCCAGCTCCCCAGTGCCACCCTGGAGGGGAGACTCACACA
 ATCTACCTTCACGCTGGAGCAACCCCTGGACCAATTCAAGAATGTCAACCTCTCTGACGCAGACCCCATC
 TGGCTGGTGGTGGCTCGCAGTGACACTGCCAGAACTTTACTGCCCCACGAAGGTGGAAGACAGGCACG
 CCCCTGCTCACTTTGACGACAATGGCTACTACCTCACACTGAGGGCCAGTCGGGTGCACTACAAGGGCAG
 CCAGAACGACAGCCGGCTTCCAGTCTCCGTGTGGGAATGACAACAATTGCTCCGTGAGATCCCAGGGC
 TGTAAATTCGCCTTTGCCAGGGCCAGGGCCCTACAGAGTGAAATTCCTGGCGATGAATGCCAAGGACCTG
 TGGCTGAGACGCTGTGGTCCAAGGACATCTACCTGCAGCAAGCCAGACATTCCGAGGAGCTCCAGGGTC
 CCAGAGCAAGGGTACTGTGGTCATTCGCTTCTCGTCAATCCTACTGGCCATCCTGCTTGCAGTCTTC
 CTCATGTTGGTCATATCTGCTTGTGCTCGAGCACTTCCGTGTCCAGCCCAGAGGAGCAAGTGCGCATGC
 GACGATATCACACGCATCATGTGGGCAGCCTCCAGCCGAGCGCAGCTCCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001109020
Insert Size:	753 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001109020.1, NP_001102490.1</u>
RefSeq Size:	975 bp
RefSeq ORF:	753 bp
Locus ID:	367994
Cytogenetics:	12q12