

Product datasheet for RN203785

Ccdc28a (NM_001037789) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ccdc28a (NM_001037789) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Ccdc28a
Synonyms:	RGD1310326
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN203785 representing NM_001037789 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAAGAGCGGAAAGCAAAGAGGAAGAGTCCCAAGTCCTTCAGTGCCCACTCTACTCAGTTGTTAATG
CCAAAAAATACCATCCCAGCTAGCAAAGCACAGGGTTTTCAAATCCCGCGTCGCAGTCAACTCCCCA
GAGGCCAAAGTTAAAAAGGGTGACGAAAGAAAAGAACAACCTCCAGGAGGAGAAGGAAGGGTGCCAG
GCCACACCCATCCAGCACTCCTTCCTCACGACGTCCTGATGTCCAGGAGATGGAGCGCGGGCTCCTTA
GTCTGCTGAATGATTTCCACTCAGGAAAGCTGCAGGCGTTTGAAATGAATGTTCCATAGAACAGATGGA
ACATGTCCGGGAATGCAGGAAAATTAGCTCGCTTGAATTTGGAGCTCTACGGGAGTTAGAAGAATT
CCGGAGGATAAGCGGAAAGCAGTCAGTGACGCCAACCTGGACAGGCTGCTGTCTGACTTAGAAGAATTGA
ATTCTTCTATACAGAACTGCATTTGGCCGATGCGCAGGACGTTCCCACTGCACCCAGCAGCTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001037789
Insert Size:	555 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001037789.1, NP_001032878.1</u>
RefSeq Size:	1267 bp
RefSeq ORF:	555 bp
Locus ID:	361454
Cytogenetics:	1p12