

Product datasheet for **RN203809**

Wrb (NM_199373) Rat Untagged Clone

Product data:

| | |
|----------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | Wrb (NM_199373) Rat Untagged Clone |
| Tag: | Tag Free |
| Symbol: | Wrb |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Cell Selection: | Neomycin |
| Fully Sequenced ORF: | >RN203809 representing NM_199373 Red=Cloning site Blue=ORF Orange=Stop codon |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGCGCGTCCGAGACAGATCGCTGGGCGTGGTTGCTGGTGCTGTGCTTCGTGTTCCGGGTGCAACGTGC
TGCATCCTTCTGCCGACCTTGTCTCGTTTATTTCCAGGGTGTGCAGAAGGATGCAGAGCAGGAGTGC
ACAGATGAGAGCAGAGATCCAGAGCATGAAGCAGGAGCTGCCACCGTCAACATGATGGATGAGTTGCC
AGATACGCCAGGCTGGAGAGAAAAATCAACAAAATGACAGACAACTCAAACCTACGTGAAGGCACGGA
CAGCTCAGTTGGCCAAGATCAAGTGTTTATAAGTGTTGCTTTCTACGTACTACAAGCAGCGTTGATGAT
CTCACTCATCTGGAAGTATTACTCTGTCCCTGTGGCTGTGGTGCCAGCAAGTGGATCACCCCACTAGAC
CGGCTGGTGGCCTTCTACTCGAGTAGCAGGTGGAATCGGAGTCACCTGCTGGATTTTAGTCTGCAACA
AAGTCGTGGCTATCATTCTCCACCCTTTCAGCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

| | |
|--------------------|--|
| Restriction Sites: | Sgfl-MluI |
| ACCN: | NM_199373 |
| Insert Size: | 525 bp |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_199373.1](#), [NP_955405.1](#)

RefSeq Size: 2419 bp

RefSeq ORF: 525 bp

Locus ID: 288233

UniProt ID: [Q6P6S5](#)

Cytogenetics: 11q11

Gene Summary: Receptor for ASNA1/TRC40-mediated insertion of tail-anchored (TA) proteins into the ER membrane.[UniProtKB/Swiss-Prot Function]