

Product datasheet for **RN203961**

Vps36 (NM_001106092) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Vps36 (NM_001106092) Rat Untagged Clone
Tag: Tag Free
Symbol: Vps36
Synonyms: RGD1309754
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN203961 representing NM_001106092
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGACCGCTTCGTGTGGACCAGCGCCTCCTGGAGATCAACGAGACCCTGGTGATCCAGCAGCGAGGAG
TGCGCGTCTACGACGGCGAGGAGAAGATAAAATTTGATGCTGGGACTCTTCTCCTTAGTACACCCGGCT
GATTTGGAGAGACCAGAAAAACAATGAGTGTGCATGGCCATTCCTATCGCAGATTGTGTTTCATCGAG
GAGCAGCGCGCTGGAATTGGGAAGAGTGCCAAAATCGTGGTTCATCTGCACCCAGCCCTCTAACAAAG
AACCCGGCCCATTCAGAGCAGTAAGAACTCTACATCAAGCTCTCCTTCAAAGAGCATGGTCAGATCGA
GTTTTACAGGCGTTTATCAGAGGAGATGACCCAAAGAAGATGGGAGACAGTGCCAGTTTTCCAGTCACTA
CAAACAAAGAAAGGACCCAGCCAGGAAGAATAAGAGCTGTTGGAATTGTAGGCATTGAAAGGAACTTG
AAGAAAAAGGAAAGAACAGACAAAAACATTTCTGAGGCCTTTGAGACCTGAGCAAGCTGATGATCAA
GGCTAAGGAAATGGTAGAGTTATCAAAATCGATCGCTAATAAGATTAAAGAGAAGCAAGGCGATGTCACA
GAAGACGAGACCATCAGGTTTAAAGTCGACTTGCTGAGCATGGGAATAGCCAACCCAGTCACCAGAGAAA
CCTATGGCTCAGGCACACAGTACCACATGCAACTGGCCAAGCAGCTGGCCGGGATATTGCAGGCTCCCTT
GGAGGAGCGTGGGGAATAATGTCACTCACAGAGGTGACTGTTTGTGTAACCGAGCCCGAGGAATGGAG
TTGCTCTCACCAGAAGACTTAGTAAATGCATGCAAGATGCTGGAAGGCCTGAAACTGCCTGTGCACTCC
GGGTTTTTCGACAGTGGTGTGATGTTAATTGAGCTTCACTCACAGGAGGAAATGGTGGCCTCGGC
CTTAGAGACGGTTTTCCGAGAGGGGATCCTTAACATCAGAAGAGTTTGCCAAGCTTGTGGGAATGTCTGTC
CTTCTGGCTAAAGAAAGGCTGCTTCTTGGCGAGAAGATGGCCACCTTTGCCGAGATGACTCTGTGGAAG
GCTTGGGTTTTATCCAAATTTATTTATGACACAGAAT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul



[View online »](#)

ACCN:	NM_001106092
Insert Size:	1161 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001106092.2, NP_001099562.2</u>
RefSeq Size:	1323 bp
RefSeq ORF:	1161 bp
Locus ID:	290851
Cytogenetics:	16q12.5