

## Product datasheet for **RN204070**

### **Dnajb1 (NM\_001108441) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Dnajb1 (NM\_001108441) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Dnajb1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN204070 representing NM\_001108441  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGCAAGGACTACTACCAGACGCTGGGCTGGCCCGGGCGCGTCCGACGACGAGATCAAGCGGGCT  
ACCGCCGCCAGGCGCTGCGCTACCACCCGGACAAGAACAAGGAGCCCGGGCCGAGGAGAAGTTCAAGGA  
GATCGCCGAGGCTACGACGTGCTCAGCGACCCGCGCAAGCGCGAGATTTTCGACCGCTATGGAGAGGAA  
GGCCTGAAGGGTGGTGGCCCAAGTGGCGGAGCAGTGGTGGTGTAAATGGTACCTTTTCAGCTACACAT  
TCCATGGAGACCCTCATGCTATGTTGCTGAGTCTTCGTTGGCAGGAATCCCTTTGATACCTTTTTTGG  
GCAGCGCAACGGGAAGAAGGCATGGACATCGATGACCCGTTCTCTAGCTTTCCGATGGGTATGGGCGGC  
TTCACCAATATGAACTTTGGACGTTCCCGCCCTACTCAAGAGCCACCCGGAAGAAGCAAGATCCCCAG  
TCACCCACGACCTTCGGGTCTCCCTTGAAGAGATCTACAGCGGCTGTACCAAGAAGATGAAAATCTCCA  
CAAGCGGCTGAACCCTGATGGAAAGAGCATTGAAATGAAGATAAGATCCTGACCATCGAAGTGAAGAGG  
GGCTGAAAGAAGGGACCAAAATCACCTTTCCAAGGAAGGGACCAGACCTCGAACAACATTCCAGCTG  
ACATCGTCTTTGTTTTAAAGGACAAGCCACACAATATCTTCAAGAGAGATGGTTCTGATGTCATTTATCC  
AGCCAGGATTAGCCTTCGGGAGGCTCTCTGTGGTTGCACTGTGAATGTTCTACTCTGGACAAGGCCTCC  
CCCTCCCCAAACACCTGAGAAACGTGGAGACCTTGTATCGAATTTGAAGTCATCTTCCCTGATAGGAT  
CCCCATCTCATCCAGAACCATCCTGGAGCAGGTTCTTCCCATATAGCTGCCTGCATTCTCAAGGACTGA  
GCCAGGGACCTTTCCAGAGCTGCAGATTTCCAGACCGTTCCACTGGTTGCATACTGCGGGAGGGCCAGG  
GAGGGCTTTGCTCTGCCGAATGTTTTCCACAGAATATATTACAATCTTTCAAAGTCGGCCTAG

**ACGGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001108441  
**Insert Size:** 1116 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001108441.1</a></u> , <u><a href="#">NP_001101911.1</a></u>
<b>RefSeq Size:</b>	2177 bp
<b>RefSeq ORF:</b>	1116 bp
<b>Locus ID:</b>	361384
<b>Cytogenetics:</b>	19q11