

Product datasheet for **RN204088**

Clba1 (NM_001008377) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Clba1 (NM_001008377) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Clba1
Synonyms:	MGC105548; RGD1307315
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN204088 representing NM_001008377 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCAAGGCGGGCGGGAGGTGGGGAGAGAGTCAGTGAGTGACCTCGCAGAGGAACTTGGAGAGGGGTCCC
 TCCACCCGACAGCCGGAGGACAGAGTGGTGATAGTTTGAACGGAGGAGGATCGGCTACGACGGTCCCGT
 GATTCTGCCGGACGCAATGCTAACAGCTCCAGGCTGGATGAGGGTCTTCTCCAGCCGTCCACACCCA
 GGGGAGCTCGGCGCGGGTGGGGAGAGTTCGAAGGCTTCAAGAGTCTTCAGCTAAGTCTGAGGAGTTTT
 CTCAGTCTTTGAACTCCTGGGAAGGGCCGCGAATGTCAGCCGCTGAGAACCCCTTCCACACCTAAAGA
 GGGCGGTTCTTGCCAGGTGCAGCAGGGTGGACCTGGGTGACAGGAACTGCTGCTGGCCCGTCTTCAGAG
 TCCATTCTCAGCTACGAGAAGGTTTTAGGCTTGCCTTCAAGAAGTCCCGGTGGAACAGGCAACCGAAG
 ACGTTTGTCTTTGGACGTTTCTTAGAAACAGGCAGCGAAGAACTGCACCTGTGCCTAGACTGTGTTT
 CGAATCTAGGAACTCTGGAGAGCCCTTCAGAACCCGACACCGTGTCTGCATCGAGGTGTCTCTGGAGT
 GAATCTCACTGCCGGGAAACCTCTTCTGTCCTTGGCATAGATGCTGCTCAGAAGAGCCTTCTGGAG
 ACCAGGTCATGATCTGGAAGTTCTGACTGCAGAAAGCCTGAAGATCTGCTGGGAGTCACTGGTTTCCA
 CCTACCACTGCAAAGCCCTGATTAGACCAAGCTCTCAGGGACATCGGGCAGCAGACAGGCGAGCCTG
 ATTACCTACAGTCTGTTTCTGAAGACCACTCCAAGGAATGGGCGCTACATCACCATCCACAGAAGA
 AGATCTTTACTCCACGCAACCTGAAAATGGCATTCTTTAATAACAACGTCTGCT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001008377


[View online »](#)

Insert Size:	966 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001008377.1</u> , <u>NP_001008378.1</u>
RefSeq Size:	1700 bp
RefSeq ORF:	966 bp
Locus ID:	362793
UniProt ID:	<u>Q5RJN9</u>
Cytogenetics:	6q32