

## Product datasheet for RN204125

### Dleu7 (NM\_001106043) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Dleu7 (NM\_001106043) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Dleu7  
**Synonyms:** RGD1305071  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >RN204125 representing NM\_001106043  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGCCAGCCCTGCGCCCTTAGTGGCTTCCATCAGCCACCAAATGGTGGCTCTGCAGACCCTGCAGCTAC  
 TACAGCAGGAGTGGGGCTGGGGCGACGGTCTAGTGCCCGAGGGAGCCCGCAGCCCGGACCACGTGCC  
 TATCGCCAGGCTCGTCGCTCAGGCCAGATGCGGACCCGCGGGGTTGGGGCGGGGAGCACGGGAGCG  
 CGGGGCTCCCCGAGGCAGGCAGCTGCGGGTGCAGAGGGGGCGCGGAGCTGCTGCCTTTTCCCCGGG  
 ACCGCGGGCCCTGTACCTTGGCGCGGATGGCGATGCGCAGCGCTTGGCTCGCGTGGTGGACTCGACCTC  
 GGAGCTGGTCAGCGTAGAGCAGACTCTGCTTGGGCCTCTCCAGCAGGAGCGGCCCTTCCCCATTCACCTG  
 AAGAGTGAGACCAGCAGCCACACACTGCTGACCTAGCTTCTTCCCATATTCATGTTCTTTGCACACA  
 GTAAGTGTGACTGCCAACTGTTCAAAGAAGACTACTCCATCCATAGTGATTAGTAAAACCAGCGCCTA  
 CAGGATAGCGTGGAAATTTAGAAACATCTGCAGTCACTTGGCTCTGCAAATTAAGGACAGCAGTTTGACA  
 GAGACT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001106043  
**Insert Size:** 639 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001106043.1</a></u> , <u><a href="#">NP_001099513.1</a></u>
<b>RefSeq Size:</b>	1693 bp
<b>RefSeq ORF:</b>	639 bp
<b>Locus ID:</b>	290308
<b>Cytogenetics:</b>	15p12