

Product datasheet for **RN204290**

Stx6 (NM_031665) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Stx6 (NM_031665) Rat Untagged Clone
Tag: Tag Free
Symbol: Stx6
Synonyms: MGC93348
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN204290 representing NM_031665
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCCATGGAGGACCCCTTCTTTGTAGTAAAGGAGAGGTACAGAAAGCAGTCAACACTGCCAGGGAT
TGTTCCAGAGATGGACGGAGCTTCTGCAGGGCCCTCCGCAGCAACGAGGGAAGAGATCGACTGGACCAC
CAACGAGCTGAGGAACAATCTCCGCAGCATAGAGTGGGATCTGGAGGACCTGGATGAAACCATCAGCATA
GTTGAAGCAAATCCTAGAAAATCAACCTGGATGCGACTGAATTGAGCATAAGAAAAGCCTTCATCACAA
GTACTCGGCAAATGTGTCAGGGACATGAAGGATCAGATGTCAGCTTCATCTGTGCAAGCGTTGGCAGAAAG
GAAAAACAGACAGGGCGCTACTAGGAGATAGCAGCAGCCAGAATTGGGATGCTGGAGTGACGGATCGATAT
GGGCGCCTTGATCGGGAGCTGCAGCTGGCCAACCTCCATTTTCATCGAGGAGCAGCAGGCACAGCAGCAGT
TGATTGTGGAACAGCAGGATGAGCAGTTGGAGCTGGTCTCTGGCAGCATCGGGGTGCTGAAGAACATGTC
CCAGCGCATCGGAGGGGAGCTGGAGGAACAGGCAGTCACTGCTGGATGATTTCTCTCAGGATGGAGAGC
ACTCAGTCTCGACTGGACAACGTGATGAAGAACTTGCAAAGTGTCTCATGACCAGTGATCGGCGCC
AGTGGTGTGCCATAGCCATCCTCTCGCGGTCTGCTGGTGTGCTGACCCTTCTCTAGTGT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_031665
Insert Size: 768 bp



[View online »](#)

| | |
|-------------------------------|---|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_031665.2 , NP_113853.1 |
| RefSeq Size: | 2515 bp |
| RefSeq ORF: | 768 bp |
| Locus ID: | 60562 |
| UniProt ID: | Q63635 |
| Cytogenetics: | 13q21 |
| Gene Summary: | involved in intracellular vesicle trafficking [RGD, Feb 2006] |