

## Product datasheet for **RN204351**

### **Gnb5 (NM\_031770) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Gnb5 (NM\_031770) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Gnb5  
**Synonyms:** MGC105255  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN204351 representing NM\_031770  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGC**C

ATGGCAACCGATGGGCTGCACGAGAACGAGACGCTGGCGTCGCTGAAGAGCGAGGCCGAGAGCCTCAAGG  
GCAAGCTAGAGGAGGAGCGGGCGAAGCTGCATGACGTGGAGCTTACCAGGTGGCCGAGCGGGTGGAGGC  
CCTGGGCCAGTTTGCATGAAGACCAGAAGGACCCTCAAAGGACACGGGAACAAAGTCTGTGCATGGAT  
TGGTGCAAAGACAAGCGGAGGATCGTGAGCTCCTCACAGGATGGGAAGGTGATTGTCTGGGATTCCTTCA  
CTACTAACAAGGAGCACGCGGTACCATGCCTTGACCTGGGTGATGGCGTGTGCTTATGCCCATCAGG  
ATGTGCCATTGCTTGTGGTGGTCTGGATAATAAGTGTTCCGTGTACCCACTGACATTTGACAAGAACGAG  
AACATGGCTGCCAAGAAGAAGTCTGTCGCTATGCACACCAACTACCTGTCTGCCTGCAGCTTACCAACT  
CAGACATGCAGATCCTGACAGCCAGCGGTGATGGCACTTGCGCCCTGTGGGATGTGGAAGTGGACAGCT  
GTTGCAGAGCTTCCATGGCCATGGGGCTGACGTCCTCTGCTTGGACCTGGCCCCCTCAGAACTGGAAAC  
ACCTTCGTGTCTGGGGTGTGACAAGAAGGCCATGGTGTGGGACATGCGCTCTGGCCAGTGTGCAAG  
CCTTTGAAACTCATGAATCTGATGTCAACAGTGTCCGATACTACCCAGTGGAGATGCTTTTGTCTCCGG  
ATCAGATGATGCCACGTGTGCCTCTATGACCTGAGGGCAGATAGAGAGGTTGCCATCTACTCCAAGGAG  
AGTATTATATTTGGAGCATCAAGCGTGGATTTCTCCCTCAGTGGCCGCTGCTCTTTGCTGGGTACAATG  
ACTATACCATCAATGTCTGGGATGTTCTCAAGGGCTCCCGGGTCTCCATCCTGTTGGACATGAAAACCG  
TGTTAGCACTCTGCGAGTGTCTCCTGATGGGACAGCATTCTGCTCAGGGTCTGGGATCACACCCTAAGA  
GTTTGGGC**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_031770



[View online »](#)

<b>Insert Size:</b>	1062 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_031770.3</a></u> , <u><a href="#">NP_113958.1</a></u>
<b>RefSeq Size:</b>	1062 bp
<b>RefSeq ORF:</b>	1062 bp
<b>Locus ID:</b>	83579
<b>UniProt ID:</b>	<u><a href="#">P62882</a></u>
<b>Cytogenetics:</b>	8q24
<b>Gene Summary:</b>	beta subunit of a heterotrimeric G-protein that is required for G-protein GTPase activity; involved as a modulator or transducer in a variety of transmembrane signalling systems [RGD, Feb 2006]