

Product datasheet for **RN204416**

Cacng5 (NM_080693) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cacng5 (NM_080693) Rat Untagged Clone
Tag: Tag Free
Symbol: Cacng5
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN204416 representing NM_080693
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGAGCACCTGTGGGAGGAAGGCCCTGACCCTGCTGAGCAGTGTCTTTGCTGTCTGTGGCCTGGCCTGC
TGGGCATCGCGGTCAGCACGGACTACTGGCTGTACTTGAAGAAGGCATAATCCTCCCCAGAACCAGAG
CACGGAGGTCAAGATGTCCCTGCACTCGGGCCTCTGGCGGTCTGCTTCCCTCGCAGGTGAGGAGCGAGGA
CGCTGTTTCACCATAGAGTATGTGATGCCCATGAACTCCCAGATGACGTCAGAATCCACAGTCAATGTTT
TAAAAATGATTCGCTCGGCCACCCCGTTCCCACTGGTCAGCCTATTCTTCATGTTTCATCGGGTTCATCCT
GAGCAACATCGGACACATCCGTCCGACCCGACCATACTGGCTTTTGTCTCTGGCATCTTCTTTATCCTA
TCCGGCCTCTCTTGTGGTGGCCTGGTGTACATCTCCAGCATCAATGACGAGATGCTCAACAGAA
CCAAGGATGCAGAGACCTATTTCAACTACAAATATGGATGGTCGTTTGCCTTTGCTGCCATCTCCTTCT
CTTAACAGAGAGTCTGGGGTGTGTCTGTGTATCTGTTTATGAAGAGGTATACTGCTGAGGACATGTAC
AGGCCCCACCCAGGCTTCTACCGCCACGGCTGAGCAACTGTTCTGATTACTCTGGCCAGTTTCTACATC
CAGATGCCTGGATCCGGGGCCGAGCCCTCGGACATCTCCAGTATGCCTCCCTACAGATGAACAGTAA
CTACCCAGCTTTGCTCAAGTGCCAGACTATGACCAGATGTCATCATCTCCCTGC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_080693
Insert Size: 828 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_080693.1 , NP_542424.1
RefSeq Size:	828 bp
RefSeq ORF:	828 bp
Locus ID:	140726
UniProt ID:	Q8VHW8
Cytogenetics:	10q32.1
Gene Summary:	member of a family of gamma subunits of voltage-dependent calcium channels [RGD, Feb 2006]