

Product datasheet for RN204426

Colec10 (NM_001130541) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Colec10 (NM_001130541) Rat Untagged Clone
Tag: Tag Free
Symbol: Colec10
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN204426 representing NM_001130541
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAATGGCTGTAGAGTCCTGCTTCGAAGAAACCAGTTCATCCTGCTGCTGCTGTTTCTCTTGCACCTTTT
 AGAGTCTGGGTCTGGATGTTGATGGTCGCTCTGCTGCAGAAGTCTGTGTACACATACCATTTCACCAGG
 ACCTAAAGGAGATAATGGTGAAGAGGTGACACAGGAGAAGAAGCAAGGATGGCAAAGTGGGACGCCAA
 GGACAAAAGGAGTTAAAGGAGAGCTGGGTGATATGGGAGCCCAGGGCGATATTGGCAAGTCTGGTCCTA
 TTGGCAAGAAGGGTGACAAAGGAGAGAAAGGTCTGCTTGGGGTACCTGGAGAAAAAGGCAAGGCAGGTAC
 CATCTGCGACTGTGGCAGATACCGAAAAGTGGTTGGACAACCTGGATATTAGTGTGCTCGTCTTAAGACA
 TCAATGAAATTCATCAAGAATGTCATAGCAGGGATCCGGGAAACTGAAGAGAAATTTACTACTATTGTGC
 AGGAGGAAAAGAACTACAGGGAATCTCTGACCCACTGCAGGATCCGGGGAGGGATGCTGGCCATGCCTAA
 GGATGAAGTTGTTAACTCTTATTGCTGACTATGTCGCAAGAGTGGCTTCTTCAGAGTGTTCATTGGC
 GTGAATGACCTTGAGAAGGAGGGCAATATGTGTTACAGACAACACTCCATTGCAGAACTACAGCAACT
 GGAAGGAGGGGAGCCTAGTGACCCTATGGCCATGAGGACTGTGTGGAGATGCTGAGCTCTGGCAGGTG
 GAACGACACAGAGTGCATCTTACCATGTATTTGTCTGTGAGTTTGTCAAGAAGAAAAAT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001130541
Insert Size: 834 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001130541.1</u> , <u>NP_001124013.1</u>
RefSeq Size:	4552 bp
RefSeq ORF:	834 bp
Locus ID:	299928
Cytogenetics:	7q32