

Product datasheet for **RN204441**

Entpd3 (NM_178106) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Entpd3 (NM_178106) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Entpd3
Synonyms:	NTPDase3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >RN204441 representing NM_178106
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTTTACGGTGTGACCCGCCAGCCCTGTGAACAGGCAGGCTTCAGGGCCCTCTCGAGGACTCCAGCCA
TCGTACACCTTGGTGGTCCTGCTTGTGAGCATCGTGGTACTTGTGACACTTACACTCATCCAGATCCACCA
CCCACAGGTTCTCTCCAGGGCTGAAGTATGGAATCGTGTGGATGCCGGCTCTCCAGAACCAGTGT
TATGTGTATCAGTGGCCAGCAGAGAAGGAGAATAACACAGGAGTGGTCAGCCAAACCTTCAGATGTAGTG
TGAAAGGCTCGGGAATCTCCAGCTATGAGAATAACCCCAAGATGCCCAAGCCTTTGAGGACTGCAT
GCTAAAGGTCAAGGAACAGGTCCCAAAGCATCTCCACGAATCCACCCGAGTTTACCTGGGGGCCACAGCT
GGGATGCGCTTGCTGAGGTTGAAAATGAGACAGCAGCTCATGAAGTCTTGAGAGCATCCAAAGCTACT
TCAAGTCCAGCCTTTTGATTTTAGGGGTGCTCAAATCATTTCTGGCAAGAGGAAGGGGTGATGGATG
GATTACAGCCAATATAATGGGAAATTCCTGGAGAAGAACCCTGGCACATGTGGGTGCACCCGCAC
GGAGTAGACACCACAGGAGCCCTGGATTTGGGCGGCGCCTCCACCCAGATATCCTTTGTGTGCGGAAGAGA
AGATGGAGCCGAACGCCAGCAGCACAGTGCAGGTGTCCTTGTACGGCTACACGTACACTCTCTACACGCA
CAGTTTCCAGTGTACGGCCGGAACGAAGCAGAGAAGAAGTTCCTGGCGATGCTTCTGCAGAGCCCTTCT
ACAGATGCCAACATCAGCAACCCTGCTACCCTCACGGTTACAGCACCACCTTCACAATGGGCCAGTGT
TTGGCAGCCTGTGCACGGAGAAGCAGAGGCCGAAGAGTACAACCCAGTGATACCATCACCTTCACAGG
AACTGGTGACCCACAGCTGTGCAGGGAGAAGGTGGCTTCTGTGTTGACTTCAGTGTGTTGCAAGAGCAA
GATGCCTGCTCCTTTGATGGCATTACCAGCCCAAGGTTCAAGGGCCATTTGTGGCATTTCAGGATTCT
ATTACACAGCCAGTCCCCTAACCTCTCTGGAAGCTTCTCCCTTACCTCCTTCAATGACAGCAGCTGGGA
CTTCTGCAGACACACTTGGAGCGAGCTCCCGTCCCTGCTCCCCAGATTTGATGAGACGTATGCCCGGTCC
TACTGCTTCTCAGCTCACTACATCTACCACTTGCTCATAAATGGATACAAATTCAGTGGCAACTTGGC
CTCAGATACGGTTTGAAGAAAGAGTGGGAATAGCAGCATCGCCTGGTCCCTAGGCTACATGCTCAGCTT
GACCAACCAGATCCCAGCTGGAAGTCCGCTGATCCAACCTACCCATACAGCCGCGGTGTTTATGGGAGTC
CTGGCCTTCTCAGGCTATCGCCTTGCTGTGCCTGGCATTCTTTTCTATCTATGTTTCAGCATTAGGA
CAAAGGAGCGCTCTGAGAATGCCTTGACCAAGCAGTGGATTCTGATTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja2034_f08.zip

Restriction Sites: SgfI-MluI

ACCN: NM_178106

Insert Size: 1590 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_178106.2](#), [NP_835207.1](#)

RefSeq Size: 3087 bp

RefSeq ORF: 1590 bp

Locus ID: 316077

Cytogenetics: 8q32