

Product datasheet for RN204447

Cenpk (NM_001106407) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cenpk (NM_001106407) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Cenpk
Synonyms:	Solt
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN204447 representing NM_001106407 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCAGAGAATCAACAACACATAGATCCAGATGTGATCGCTGCTAAAGATACAGAGGAAGAAGCTTATCA
AACAAATGCGAAGAAATGTGAAAGACGTGGAAGAATGTCAGAATAAATTACTTATTGGACCCGAGAC
ACTCACTAACTCAATGCTCAGCTCTCATTACTTATTATGCAAGTGAAATGTTTGACTGCTGAAGTTGGT
CAATGGAACAAAAGAAAACCCGAAACAGTCCCCTCACTGAAGGCGTTTTGTTAACATTAGGAAAAGAAG
AGTTCCAAAAGTTGAGATATGATCTTGAATGGTACGGTCTATCGTTCAATCAGAGAATGAAAAGTTAAA
AGAAGACTTAGCAAGGGAGCAACAGTGGTTGGATGAACAGCAGCAAACTGGAAGCTTAAATGTATTA
CACAGTGACTTGAATAATCGGGTTGTAACATTTACTGAATCAAGAATTTTAAATGAAGTAAAGCCAAGA
TCCTTAGTATAAAAAGAAATTTAAAGAGAAGCTCCTACTTGCCTGGGTGATTTTTTGAAGATCATTTTCC
ACTGCCTGATGGAATACTCAAAAGAAAAGGAAAAATTTCAAGAATTGAATACACAGCTGATAACACTG
AATGAAATGTTAGAGGTCCCACTCCTGGACATGGGATCTACCATTGAGAGTTGTGTGTTCCAGTGAGA
CTTTCTTGAAGAAAACCAAACTTCATTTGCAAGTGCTGCAATTGGAGACAGCTTCTGGATTAGGGGTGG
GTTCCCTCCTAGTTCTAGGACTAGTAAATCTAGAAGAGTATTTTCAGAACTCAATGAGTTGGTTATAAGC
AATCAGTTTGTACCTAAGAATCTCTACTAATAAATTTTTTATATAAATATCTACTGTTTGA

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001106407
Insert Size:	909 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001106407.1, NP_001099877.1</u>
RefSeq Size:	1641 bp
RefSeq ORF:	909 bp
Locus ID:	294712
Cytogenetics:	2q12