

Product datasheet for RN204451

Pthlh (NM_012636) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pthlh (NM_012636) Rat Untagged Clone
Tag: Tag Free
Symbol: Pthlh
Synonyms: PLP
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN204451 representing NM_012636
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGCTGCGGAGGCTGGTTCAGCAGTGGAGCGTCCTGGTGTTCCTGCTCAGCTACTCCGTGCCCTCCCGG
 GCGCTCGGTGGAGGGCTTGGTCGCAGGCTAAAACGCGCGGTGTCTGAGCACCAGCTACTGCATGACAA
 GGGCAAGTCCATCCAAGACTTGCGCCGCGTTTCTTCTCCACCATCTGATTGCGGAGATCCACACAGCT
 GAAATCAGAGCTACCTCGAGGTGTCCCTAACTCCAAACCTGCTCCCAACACAAAAACCCCTGTGC
 GTTTTGGTTCAGACGACGAGGGCAGATACCTAACTCAGGAAACCAACAAGGTGGAGACGTACAAAGAGCA
 GCCACTCAAGACGCCCGGAAGAAGAAGAAAGGCAAGCCGGGAAACGCAGAGAACAGGAGAAAAAAG
 CGAAGGACTCGGTCTGCCTGGCCAGGCACAAGTGGGAGTGGCCTGCTTGAGGACCCAGCCACACCT
 CCCCAGCTCGACCTCCCTGGAGCCAGCTCAAGGACGCAT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_012636
Insert Size: 534 bp



[View online »](#)

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_012636.1 , NP_036768.1
RefSeq Size:	1087 bp
RefSeq ORF:	534 bp
Locus ID:	24695
UniProt ID:	P13085
Cytogenetics:	4q44
Gene Summary:	induces vasodilation; may play a neuroprotective role in response to ischemic brain injury [RGD, Feb 2006]