

## Product datasheet for **RN204452**

### Suv39h1 (NM\_001106956) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Suv39h1 (NM\_001106956) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Suv39h1  
**Synonyms:** RGD1565028  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN204452 representing NM\_001106956  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGGGAACCGCGACTCTAGGTTGCAGTGTGTGCTGTAATCTTCTGGAATCAGCTGCAGGATCTGT  
GCCGACTAGCCAAGCTTTCTTGTCTGCCCTTGGTGTCTTAAGAAGAATCTGTATGACTTTGAAGTTGA  
ATACCTGTGTGATTATAAGAAGATCCGTGAGCAGGAGTATTACCTGGTTAAATGGCGTGGGTATCCTGAC  
TCAGAGAACACCTGGGAGCCACGGCAGAATCTCAAATGTGTGCGATTCTTAAGCAGTCCACAAGGACT  
TAGAAAAGAGAGCTTGTCCGGCGACACCGCCGGTCAAAGCCACCCAGGCATCTGGACCCAAACTTAGCCAA  
TTACCTGGTGCAGAAGGCCAAGCAGAGGCGGGCACTGCAGCGTTGGGAACAAGAGCTCAATGCCAAGCGC  
AGCCATCTGGGGCGGATCACTGTGGAGAATGAGGTAGACCTGGATGGCCCTCCAAGGTCCTTTGTCTATA  
TCAATGAGTATCGAGTTGGTGGGGCATCACCCCAACCAGGTAGCTGTTGGCTGTGAGTGCCAGGACTG  
TCTGTTGGCACCCACTGGAGGCTGTTGCCCTGGAGCATCCCTGCACAAGTTTGCCTACAATGACCAAGGC  
CAGGTGCGACTGAAAGCTGGGCAGCCCATCTACGAGTGCAACTCTCGCTGTTGCTGTGGCTATGACTGCC  
CAAACCGTGTAGTCCAGAAAGGCATCCGCTACAACCTCTGCATCTCCGCACTGATGATGGCCGAGGCTG  
GGGTGTCCGACGCTGGAAAAGATCCGCAAAAATAGCTTTGTTATGGAGTATGTGGGAGAGATTATTACC  
TCAGAGGAGGCAGAGCGGAGGGCCAGATCTACGACCGCCAGGGCGCCACCTACCTCTTTGACCTGGACT  
ACGTGGAGGACGTATATACCGTGGATGCCGCCTATTATGGCAACATCTCTCACTTTGTCAACCATAGTTG  
TGATCCCAACCTGCAGGTGTACAACGTATTATAGACAACCTTGATGAGCGACTACCCCGCATCGCATTC  
TTTGCCACAAGAACCATCTGGGCGGGCGAGGAGCTCACCTTTGATTACAACATGCAAGTGGACCCCGTGG  
ACATGGAGAGTACCCGAATGGACTCCAACCTTTGGCCTGGCAGGGCTCCCCGGCTCCCCCAAGAACGGGT  
CCGTATTGAATGCAAATGTGGGACAACAGCTTGCCGAAAATACCTCTTCTAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online >](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_001106956
<b>Insert Size:</b>	1242 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001106956.1</a></u> , <u><a href="#">NP_001100426.1</a></u>
<b>RefSeq Size:</b>	3337 bp
<b>RefSeq ORF:</b>	1242 bp
<b>Locus ID:</b>	302553
<b>Cytogenetics:</b>	Xq12