

Product datasheet for **RN204482**

Aqp11 (NM_173105) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Aqp11 (NM_173105) Rat Untagged Clone
Tag: Tag Free
Symbol: Aqp11
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN204482 representing NM_173105
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCCGCGCTACTGGGGCTCCCGCCGAGGTGCAGGACACCTGCATCTCGCTGGGACTAATGCTGTTGG
TCGTGCTGTTTCATGGGGCTGGCCCGCGTGATCGCCCGCAACAGCTGCACAGGCCCATGGTCCACGCCTT
CGTCCTGGAGTTTCTAGCTACCTCCAGCTCTGCTACTGCACCCATGAGCTCCAAGTCTGAGCGAGCAG
GACTCGGGGCACCCACCTGGACGCTGACTGATCTACTTCTTTTCTTGGTGCATGGCCTGACCCTGG
TGGCACAGCTAGCAACCCGTGCGCGTGATGATGCAGATGATTCTGGGGGCATGTCCCCGAGATGGG
TGCGGTGAGGTTGATGGCTCAGCTGTTAGTGCCCTATGCAGCAGGTAAGTGCATTAAGTGCCTGTGGAGC
CTGAGCCTGACCAAGTACCATTTGACAAGAGGATTTAGCTTGCAGGAATCCCATCAACACCGACATAT
CCAAAGCGATCATCATAGAAGCCATCTGCTCCTTTATTTTCCACAGCGCTCTACTGCACTCCAGGAGGT
CCGAACCAAGCTTCGCATCCACGTGCTGGCTGCACTCATCACCTTTTGGCCTATGCAGGAGGGAGCCTC
ACAGGAGCGCTTTAAACCAGCACTGGCACTTTCGCTGCACTTCCCGTCTTTGACGAATCGTTCTATA
AGTTTTTCGTAGTATATTGGGTTGCTCCTTCTAGGCGTACTGCTGATGATCCTCATGTTCAAGCTTCT
CCTCCGTGGCTGCACAACAACCAACTGAGTAACAAGAAGGA**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja3129_e03.zip
Restriction Sites: SgfI-MluI
ACCN: NM_173105
Insert Size: 816 bp



OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_173105.1](#), [NP_775128.1](#)

RefSeq Size: 1299 bp

RefSeq ORF: 816 bp

Locus ID: 286758

UniProt ID: [Q8CHM1](#)

Cytogenetics: 1q32

Gene Summary: member of the major intrinsic protein family [RGD, Feb 2006]