

## Product datasheet for **RN204509**

### Adamdec1 (NM\_001106046) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Adamdec1 (NM\_001106046) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Adamdec1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN204509 representing NM\_001106046  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGCCTCTACAAAGCACAGATGAGGGGAACATGCTGTTATTCCACACAACCTGGGAAGAACCAGGTGCAG  
TTAACTACAAGCACAGTGAGAAGCAAATTGGCAGGAAACGAAGCCATCTCGAACTCCAGGTCACTCAA  
AAGCCAAATGTAACACTCACTTCAAGGACAGAAGTACATCGGTCTCTTTTGGTGTGGATAATTCCTAT  
TATAAGATGTATAAAGGGAATGTGACTCGAATGAGAAGCTTTGTGTTGAGGTGCTGAATCTGCTCAACG  
TGATTTATAAAACCATAGATATTCAAGTGTCTTTGGTGGCATGGAGATCTGGTCAGACCCGGATAAGAT  
AAAGGTGGTGCCCAACCTTGGTGCCACATTACAAAATTCATGAGATGGCATTATTCTAACCTGGGACAA  
AGGATCCATAACCATGCCAGCTTCTCAGTGGAGCTGGCTTCTCCATGGACGTGTAGGAATGCCAGCTG  
CGAATTCCTGTGACTACTACATCTCAGTGTCTGTAATTGAGGCTAAAAGAAAGAATAATGTGGCTTTGT  
AGCACTGATGTCACATGAGTTGGCCATGCCCTTGGTATGCAGGATGTTCCATATTATACCAAGTGTCC  
TCTGGCAGCTGTGTGATGAATCAGTACCTGAGCTCAAAATCCCCAAGGATTTCACTACTGTCAGCCGCT  
CACGGTTTCAAGGATTCCTTTCATCTAGAAATACAAGGTGTCTGCTGTTGGCTCCTGACCCTAAAAACAT  
AATAAAACCAACATGTGGAACCGACTTCTGGATATGGGAGAAGGCTGTGACTGTGGCTCTCCTGAGGAA  
TGTACCAATCTTTGCTGTGAGCCCTGACATGTAGGCTGAAATCTCAACCTGACTGCAGAGAGGAAGCTT  
CCAACCTCGTCATAGAA**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001106046  
**Insert Size:** 930 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001106046.1</a></u> , <u><a href="#">NP_001099516.1</a></u>
<b>RefSeq Size:</b>	2394 bp
<b>RefSeq ORF:</b>	930 bp
<b>Locus ID:</b>	290338
<b>Cytogenetics:</b>	15p11