

Product datasheet for **RN204697**

Capns2 (NM_001109380) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Capns2 (NM_001109380) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Capns2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN204697 representing NM_001109380 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTTCTCGCAAAGGCCATCCTGGAAGGTCAGATCGGGGTCTCGGAGGGGCTCTTGAGGCCTTCTTG
 GAGGAGGTGGTCAGGCAGGAGCAGGTAGGGGCAACATCGGAGGGATCATCGGAGGAATTGTCAATTTTAT
 CAGTGAGGCCGCCGAGCTCAGTACACCCAGAGCCGCTCCCCACAGCGGCATTTCACTGTTGTGGAG
 GCCTCAGAAAGTGACGAAGTCAGGCGTTTCCGGCAGCAATTTACACAGCTGGCTGGACCCGACATGGAGG
 TGGGTGCCGCCGACCTGATGAACATTCTCAACAAGGTCCTCTCTAAGCACAAGGAGCTGAAGACGGACGG
 TTTCAGCCTTGACACCTGCAGGAGCATCGTGTCTGTCTCATGGACAGTGACACGACTGGGAAACTGGGCTTT
 GAGGAATTTAAGTACCTGTGGAACAACATCAAGAAATGGCAGTGTGTTTTTAAGCAGTATGACAGCGACC
 ATTCTGGCTTCTCAGGAGCTCCAGCTACACGGGGCCATGCAGGCAGCAGGCTTCCAGCTAAATGAGCA
 GCTTTACCTAATGATTGTCCGTCGGTATGCTGACGAAGACGGCAGCATGGACTTTAACTTTATCAGC
 TGCCTGGTTCGCCTGGATGCCATGTTTCGAGCTTTCAAGACTCTGGATCGTGATAGGGATGGGCTGATTC
 GGGTTTCCATCAGGGAATGGCTGCAGCTGACCATGTATTCT**GATGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001109380
Insert Size:	744 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001109380.1, NP_001102850.1</u>
RefSeq Size:	943 bp
RefSeq ORF:	744 bp
Locus ID:	679870
Cytogenetics:	19p11