

Product datasheet for **RN204774**

Pabpc6 (NM_001106208) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pabpc6 (NM_001106208) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Pabpc6
Synonyms:	Pabpc3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF:	>RN204774 representing NM_001106208 Red=Cloning site Blue=ORF Orange=Stop codon TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCCCGGATCGCC ATGAACCCAGCGACCCAGCTACTCCTTGGCTTCACTCTATGTGGGGACCTGCACCCCGACGTGACCG AGGCGATGCTCTACGAGAAATTCAGCCCGCAGGGCCCATCCTCTCCATCCGGGTGTACAGGGACAGAAC CACCCGCCGCTCGCTGGGCTACGCGTCTGTAACTTCCAGCAGCTGGAGGATGCCGAGCGGGCCTTGAC ACCATGAATTTGATGTGATAAAGGGCAAGCCGGTTCGCATCATGTGGTCTCAGCGGGACCCATCGCTGC GCAAGAGTGGGTAGGCAACATATTCGTCAAAAATCTGGACAGATCCATTGACAGTAAACGCTGTATGA TACGTTTTCTGCCTTCGGTAACATCCTTCTGTAAAGTGGTGTGTGATGAGAATGGTTCCAAGGTTAT GGATTCGTACACTTCGAGACACAGGAAGAAGCCGAACGAGCAATTGAGAAAATGAACGGGATGTTCTTAA ATGACCGCAAAGTGTTCTGGGCGCTTAAAGTCTCGTAGAGACCGGCAGGCTGAACTTGGAGCCAGGGC CAAGGAGTTCACCAACGTTTACATCAAGAATTTGGGAGAAGACATGGATGATGAAAGACTCCAGGATCTC TTCAGCAGATTTGGGCTGCCTTAAAGTGTAAAGTCATGACTGATGAAAGTGAAAGTCCAAAGGATTCCG GATTTGTAAGCTTTGAAAGGCATGAAGATGCACGAAAGGCAGTGGATGAAATGAACGGAAAAGACCTCAA TGGAAAACAAATTTATGTAGGTCGAGCTCAGAAAAAGGTGGAACGACAGACTGAACTGAAGCACAAATTC GGGCGATGAAGCAAGATAAGCCCAAAATTAACAGGTACCACAAGATAGAAGCGTCAGATGCCAAGGTG TTAACCTTTACGTGAAAACCTCGATGATGGTATTGATGACGAGCGTCTCCGGAAGGAGTTTTCTCCGTT TGGGACTATCACTAGTGCAAAGGTAACGATGGAGGGTGGGCGCAGCAAAGGATTTGGTTTTGTATGTTTC TCCTCTCCAGAAGAAGCCACCAAGCCGTTACAGAAATGAATGGCAGAATAGTGGCAACCAAACTGT ACGTAGCTTTAGCTCAACGCAAAGAGGAGCGCCAGGCTCACCTCAGTAACCAAGTACATGCAGAGAATGGC AAGCACCAGTGCCGTGCCCAACCCAGGCATCAACCCTTTCCAGCCAGCACAGGCTCTTTCAGGATACTTC ATGACAGCTACAGCACCGACTCAGAACCGTGGGACATACTATGCTCCTAACCAAACTACTCAGCAAGGAC CAAGTGCCCGTGGGAGTGCTCAGGGTACCAGAGCGCATCCATTCCAAACATGTCTAGTGCCATCCACCC ATCTAATACTATGCCATCATTTAGTACTTCAGGACCCACTTCTCAGAGCTGTACGACACCGTACTGCC AGCACATCCACACAGATGATGGGTCCACATCCTGCAGTGCTTCTACTGCAGTGCTCCTGCCGCCCGCA CAGTTACTCAGTATAATATACTGCTGGAGTCCGAAATCCCCCTCAGCATCTTAATACACAGCCACAGGT CACCACGACGCGTCTGCTGTTTCTGCCAAGGTAAAGAATCTTGACTGTCTCCATGTTGGCATCTGCT CCTCCTCAAGCTCAAAAGCAAAATGTTAGGGGAATGGCTATTTTCTTATTACAGGCCATGCACCTGCC TTGCCGTTAAATCACCGGTATGCTGTTGGAAATTGATAATTTAGAGTACGCCACATGCTTGAGTCTCC AGAGTGTCTCCACACCAAGTCGATGAAGCCATAGCTGTACTGCAAGCCACCAAGCTAAAGAGACTTCC CAGAAAGCAGTCAGCAGTGTTGCTGGTGTCCCACTTCTTAA AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC TGGATTACAAGGATGACGACGATAAGGTTTAA
Restriction Sites:	Sgfl-RsrII
ACCN:	NM_001106208
Insert Size:	1932 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001106208.2, NP_001099678.1

RefSeq Size: 3100 bp

RefSeq ORF: 1932 bp

Locus ID: 292295

Cytogenetics: 1q12