

Product datasheet for **RN204830**

Tgfa (NM_012671) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tgfa (NM_012671) Rat Untagged Clone
Tag: Tag Free
Symbol: Tgfa
Synonyms: RATTGFAA; TGFAA
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN204830 representing NM_012671
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGTCCCCGGCCGGACAGCTCGCTCTGCTAGCGCTGGGTATCCTGGTAGCTGTGTGTCAGGCTCTGG
AGAACAGCACGTCCCCCTGAGTGACTACCCGTGGCGGCTGCAGTGGTGTCTCACTTCAACAAGTGCC
AGATCCCACACTCAGTATTGTTCCATGGGACCTGCAGGTTTTTGGTGCAGGAAGAGAAGCCAGCATGT
GTCTGCCACTCTGGGTACGTGGGTGTTGCTGTGAGCATGCAGATCTCCTGGCAGTGGTGGCTGCCAGCC
AGAAGAAGCAGGCCATCACTGCCCTGGTGGTAGTCTCCATTGTGGCCCTGGCTGTCTCATTATCACCTG
TGTGCTGATCCACTGCTGTGAGTCCGCAAGCACTGTGAGTGGTGGCCGTGCCCTCGTCTGCAGACACGAG
AAGCCCAGCGCCCTCTGAAAGGAAGGACTGCCTGCTGCCACTCTGAGACAGTGGCT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_012671

Insert Size: 480 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_012671.2, NP_036803.1</u>
RefSeq Size:	4222 bp
RefSeq ORF:	480 bp
Locus ID:	24827
UniProt ID:	<u>P01134</u>
Cytogenetics:	4q34
Gene Summary:	acts as an EGF receptor ligand; inhibits apoptosis that results from epithelial detachment from the extracellular matrix [RGD, Feb 2006]