

Product datasheet for **RN204877**

Tifab (NM_001025029) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tifab (NM_001025029) Rat Untagged Clone
Tag: Tag Free
Symbol: Tifab
Synonyms: RGD1311681
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN204877 representing NM_001025029
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAGAGGCCCTCACAGTCTGCAGGTGAGCCTGTACCACCCGACACAGGGCCCGGTGCGCTTTGCCA
AAGTTCGGCTGCAGCTGCAACACGACGCCAGCCGGTTGCTGGTTGGGCGAGGGCAGGACACCCACCTCCA
GCTGCAGCTGCCACAGTTGTCCCGCCACCATTTGTGCCTGGAGCCCTACCTGGAGAAAGGCAGCAACCTG
CTGGCCTTCTGCCTCAAGGTGCTGACCCGCCAGAGCTGTGTGGGTCAACGGTTGCCGCTGAGGTACC
TGGAGCAGGTTCCCTGCACTCCGTCAACAGAATCTCCTTCTCTGGCATCCAGATGCTCATCCGCAAGA
GGGGGGCGCCTCCCTTGAGGCCTTTGTCTGCTATTTCCATCTCAGCCCTCGCCCTGATTTACAGACCC
AAGGCTCAGGAGACAGACGAATGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001025029
Insert Size: 444 bp
OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001025029.1</u> , <u>NP_001020200.1</u>
RefSeq Size:	1377 bp
RefSeq ORF:	444 bp
Locus ID:	364674
UniProt ID:	<u>Q5BK67</u>
Cytogenetics:	17p14
Gene Summary:	Inhibits TIFA-mediated TRAF6 activation possibly by inducing a conformational change in TIFA. [UniProtKB/Swiss-Prot Function]