

## Product datasheet for **RN204889**

### Dap3 (NM\_001011950) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Dap3 (NM\_001011950) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Dap3  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN204889 representing NM\_001011950  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGAGAACAAGACCATGCTGACAGGAATGACAAGTTTTTCTCCAGGGTCCAGAAGTTGGACCCGAGAT  
GTTTTTGCACGTGAGTGTCCAAGCACCTCAGAACAATCAGGTTCCAGCTGAGCGTCCGAGAGCTGTTTC  
CCGCACCACTGACAGTGACCCAGCTAAGCACGGGAGCAGCATGAGGGTCAGCACTACAGCATCCCCCA  
CAGGATCTGAAGACGGTGTTCCTCACGGCCTGCCTCCTCGATACGCAATGCAGGTGAAGACTTTTGGTG  
AAGCTTGCTGATGGTTAGGAAGCCAGCTCTCGAGCTTCTGGGCTACCTGAAAAATACCAATTTTGTCTA  
TCCAGCTGTACGATATCTTCTATGGAGAGAAGGGGACAGGAAAGACCCTCAGTCTCTGTCTGCTGTT  
CATTTCTGTGCAAAAAGTACTGGCTGATTCTGCACATCCAGATGCTCATCTTTGGGTAaaaaactGTG  
GGGAGCTGCTACAGTCCACTTACAACAACAGCGCTTTGACCAACCCTTAGAGGCCTCTACCTGGCTGAA  
GAATTTCAAACTACCAACGAGCGCTTTCTGGGTGAGTAAAAGTTCAAGAGAAGTATGTTTGAATAAG  
AGAGAAAGCACTGAGAAAGGCAGTCTCTGGGAGAAGTTGTGCAACAGGGCCTAACGCGAGTGAGAAACG  
CCACTGACGCTGTTGGGTTGTGCTGAAGGAGCTGAAGACGCAGAGCTCCCTAGGACTCTTTCACCTCT  
GGTGGCTGTGGATGGAGTCAACGCCCTCTGGGGAAGGACCACGCTGAAAAGAGAAGATAGAACCTTGATT  
GCCCCGGAGGAACTCGCCCTGTCCGCAACCTAAGAAAAATGGTAAAAATGACTGGCACGGAGGAGCAA  
TTGTGCTGAGCTTGAACCAACCGGGTCTCTCTCAAGTCCAGAACAGCCTATTTGCCACATGAGCTGCT  
GGGAAAGGAAGGATTTAATGCCCTGGAGCCTTTTCTTCCATCCTCATTCCAACATAATCCAAAGGAA  
TTTGAAAGTTCCTTTCAGTATTACTTGGAGAACAAGTGGCTTTCAGCATGAGAAAGCTTCAACAGAGGAA  
GGAGGAAGGAGCTGAAGTTCCTGAGTAACTGCAACCCTGAGCAGCTAGAGCGGCTCTGCGCTCACTGTG  
A

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI



[View online »](#)

ACCN:	NM_001011950
Insert Size:	1191 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
RefSeq:	<u>NM_001011950.2</u> , <u>NP_001011950.2</u>
RefSeq Size:	1870 bp
RefSeq ORF:	1191 bp
Locus ID:	295238
Cytogenetics:	2q34