

Product datasheet for RN204906

Igfbp3 (NM_012588) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Igfbp3 (NM_012588) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Igfbp3
Synonyms:	IGF-BP3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN204906 representing NM_012588 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCATCCCGCGCGCCCCGCGCTCTGGCGGGTGCCTCACCGCCCTCACTCTGCTCCGCGGACCGCCAG
TGGCGGGGCGCGGGCGGGTGGCGCGGGCCCCGTGGTGCCTGCGAACCGTGCGACGCGCGTGC
GCTGGCCAGTGCAGCCTCCGCCACCGCGCCCGTGCACGGAGCTGGTGCAGAACCCGGCTGCGGC
TGCTGCCTGACTGCGCGCTGCGGAAGGCGACGCTGCGGCGTCTACACGGAGCGCTGTGGACCGGCC
TCCGCTGCCAGCGCGACCGGCCGAGCAGTATCCCTGAAGGCGCTGCTGAATGGCCGCGGTTCTGCGC
CAACGCCAGCGCCGCGCAGCAGCCTGAGTGCCTACCTCCCCTCCCAGCCGTCTCCTGGAAACACCACTGAG
TCTGAGGAGGACCACAATGCTGGGAGTGTGAAAGCCAGGTTGTCCCCAGCACACACCGCGTACTGATT
CCAAGTTCATCCACTCCATTCAAAGATGGAGGTCATCATAAAAGGCCAGGCTAGGGACAGCCAGCGCTA
CAAAGTTGACTATGAGTCCCAGAGCACAGACACCCAGAATTCTCCTCCGAGTCTAAGCGGGAGACAGAA
TATGGTCCCTGCCGAGAGAAATGGAGGACACTGAATCATCTGAAGTTCCTCAATGTGCTGAGTCCCA
GGGCGTCCACATCCCAAAGTGTGACAAGAAGGGGTTCTATAAGAAGAAACAGTGTGCCCTTCAAAGG
CAGAAAGCGGGGCTTCTGCTGGTGCCTGGACAAGTACGGGACGCCATTGCCAGGCTATGACACCAAGGGG
AAAGACGACGTGCATTGCCTCAGCGTGCAGAGCCAG**TAG**

AC**GGGCGCCG**CTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGA
TTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-NotI
ACCN:	NM_012588
Insert Size:	879 bp



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_012588.2](#), [NP_036720.2](#)

RefSeq Size: 2416 bp

RefSeq ORF: 879 bp

Locus ID: 24484

UniProt ID: [P15473](#)

Cytogenetics: 14q21

Gene Summary: an insulin growth factor (Igf) binding protein; involved in modulating IGF-I action [RGD, Feb 2006]