

Product datasheet for **RN204913**

Car13 (NM_001134993) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Car13 (NM_001134993) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Car13
Synonyms:	RGD1560453
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN204913 representing NM_001134993 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGCGAGGTTAAGCTGGGGATACGACGAGCACAAACGGTCCTATTCAGTGAATGAATTGTTCCCTATTG
CGGATGGTGATCAGCAGTCTCAATTGAGATTAACCAAGAAGTAAATACGACTCCTCACTCCGACC
TCTCAGTATCAAGTATGATCCTGCCTCGGCTAAATCATCAGCAACAGTGGCCATTCTTCAACGTTGAC
TTTGACGACACAGAGGACAAATCAGTTCTGCGTGGAGGTCCTCTCACTGGAAGCTACAGGTTGCGGCAGT
TCCATCTGCACTGGGGTACGACGACACCATGGCTCTGAGCATGTGGTGGATGGCGTGAGGTATGCTGC
AGAGCTGCATGTTGTCCACTGGAATTCAGACAAATACCCAGTTTTGTAGAGGCAGCTCATGAGTCCGAT
GGGCTGGCTGTCCTGGGAGTATTTCTACAGATTGGGGAACACAATCCTCAACTGCAAAAGATCACGGACA
TTTTGGATTCCATTAAGGAAAAGGGTAAACAAACACGGTTCACAAATTTGACCCCTTATGTTTGCTTCC
GTCATCCTGGGACTACTGGACATACCCTGGTCTCTGACGGTTCACCTCTTCTTGAGAGTGTACAGTGG
ATTGTTTTAAACAACCTATAAGTATCAGCTCTCAACAGTTGGCCAGATCCGAAGCCTCTTGTCACAG
CTGAGGGCGAGAGTGGGCTTTCTGTTAAGCAATCACCGTCCGCCACAGCCCTGAAGGGCCGACAGT
GAGAGCATCTTTCTATTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001134993
Insert Size:	789 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001134993.1, NP_001128465.1</u>
RefSeq Size:	1982 bp
RefSeq ORF:	789 bp
Locus ID:	499566
Cytogenetics:	2q23