

Product datasheet for **RN204953**

B4galt7 (NM_001031661) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: B4galt7 (NM_001031661) Rat Untagged Clone
Tag: Tag Free
Symbol: B4galt7
Synonyms: MGC116242
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN204953 representing NM_001031661
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTTGCCCTCCCGAGGAAAGCGGCGCAGCTGCCCTGGGAGGATGGCAGGGCCAGGCTGCTTCCTGGAG
 GCCTCCGCGGAAATGCTCCATCTCCACCTCTTCATTGCCTTTCTCCTGTTGGTCTTCTCTCCCTGCT
 CTGGCTGCAGCTCAGCTGTTCTGGAGATGTGGCCAGGTGACCAGGGGACAAGGGCAGGAGACCTCGGGT
 CCACCCGAGCTTGTCCCCAGAACCGCCCTGAGCACTGGGAAGAAGATGAGTCATGGGGCCTCACC
 GGTGGCAGTGTGGTGCCTTTTCGGAACGCTTTGAGGAGCTGCTGGTCTTTGTGCCCCACATGCACCG
 CTTCTAAGCAGGAAGAAGATCCAGCACCATCTATGTGCTCAACCAGGTGGATCACTTCAGGTTCAAT
 CGGGCAGCACTCATCAATGTGGGCTTCTGGAGAGCAGCAACAGCACAGACTACATCGCCATGCACGATG
 TGGACCTACTCCCTCTCAATGAGGAGCTGGACTATGGCTTCCCGAGGCTGGGCCCTTCCATGTGGCTTC
 CCCAGAGCTCCACCCTCTACCACTACAAGACCTATGTGGTGGCATTCTGCTGCTGTCCAAACAGCAC
 TACCAGCTGTGCAACGGGATGTCAACCGCTTCTGGGGCTGGGGCCGAGAGGACGATGAATTCTACCGGC
 GCATCAAAGGAGCTGGTCTCCAGCTTTTCGCCCCTTGGGGATCACAACTGGGTACCAGACATTTGCCA
 CTTGCATGACCCAGCCTGGCGGAAGAGGGACCAGAAACGATTGCGGCTCAAAAACAGGAACAGTTCAAG
 GTGGACCGGGAGGGAGGCTGAACACTGTGAGGTACCGGGTGGATTCCCGCACGGCACTGTCCGTAGGAG
 GGGCCCCCTGCACTGTTCTCAATGTCATGCTGGACTGCGATAAGACAGCCACTCCATGGTGCACATTTGG
 CTGA

ACGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001031661



[View online »](#)

Insert Size:	984 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001031661.1</u> , <u>NP_001026831.1</u>
RefSeq Size:	1591 bp
RefSeq ORF:	984 bp
Locus ID:	364675
Cytogenetics:	17p14