

Product datasheet for **RN204966**

Mak16 (NM_001014002) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mak16 (NM_001014002) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Mak16
Synonyms:	Rbm13; RGD1311297
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN204966 representing NM_001014002 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGCAGTCGGACGATGTTATCTGGAATACACTAGGAAACAAACAATTTTGTCTTACAAAATAAGAACCA
 AGACTCAGGGCTTTTGCAGAAATGAATATAGTCTGACAGGTCTGTGCAACCGGTCGTCTGTCCCCTGGC
 AAATAGTCAGTATGCTACCATTAAGAAGAGAAAGGGCAGTGCTACCTGTACATGAAGGTCATAGAGCGC
 GCAGCCTTCCCTCGCGTCTCTGGGAGCGGGTCCGGCTTAATAAAAACTATGAGAAAGCACTGGAACAAA
 TAGATGAAAATTTGATTTACTGGCCCCGCTTTATCCGACACAAATGTAAGCAGAGGTTTACCAAGATAAC
 CCAGTACCTGATTCGGATTAGAAAACCTTACACTGAAGCGAAGGAAGAACTTGTTCTTTGAGTAAGAAG
 GTTGAAAGAAGGAAAAACGAAGAGAGGAAAAAGGCATTATAGCCGCCAGCTGGACAATGCCATTGAGA
 AGGAATTGCTGGAGAGACTGAAACAAGACACGTATGGCGACATCTACAACCTCCCATCCACGCCTTCGA
 CAAGGCCCTGGAGAAACAGGAAGCAGAAAGTACTCTGAAGATGAAGAAGATGAGGAGGAAGATGAAGAA
 GAGGATGTGGGAAAAGAGATTTGTAGAAGATGATGAGGTAGAGGAAAGTACCTGAGTGACTTTGAGG
 ATATGGATAAGCTGAACACTGACAGTGAGGAAGATGATGAATCCTCCCGTGAAGAAGAGGCTCATAAGGC
 CAAACACAAAGGCAAAGCACCTTGAAAGGACCTTTGCGGAAAAAACGGGCTATGTGGAATAGAGTAT
 GAACAGGAGACAGAGCCCATGGCCAAAGCCAAAGCCACA**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001014002
Insert Size:	882 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001014002.1, NP_001014024.1</u>
RefSeq Size:	1741 bp
RefSeq ORF:	882 bp
Locus ID:	306526
Cytogenetics:	16q12.3