

Product datasheet for **RN204980**

Rangrf (NM_001105790) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rangrf (NM_001105790) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Rangrf
Synonyms:	RGD1563195
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN204980 representing NM_001105790 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAGCCCGCCAGAACTGTCCACTGTTTCGGAGGCGCCTTTTCCGCCATCCTCCCCGAGGGGCCATTG
 ATGTGAGTGACCTCCGACCGGTGCCAGACAACCAAGAAGTTTCTGCCATCCTGTGACCGACAGAGTCT
 GATCATAGAATTCTGGAGCTGCAGGCCACGTGCAGGGCGAAGCGGCTGCGCGGTACCATTTTGAGGAC
 GTTGGCGGGGTGCAGGGGCTAGGAGTGTGCACGTGCTATCTGTGCAGCCTCTCTGTTTGGAGAATTAT
 CCCTGAGAGGCTGCTGTCAAGATGCCTGGTTCTTTCTGGCAAGCAGCAGGTAGCTAAAGAAAACCAACA
 GGTAGCAAAGGATGTCACTGCATCAGGCCTTGCTTCGGTTGCCCAAGTATCAGACTGATCTTTTGCTC
 ACCTTCAATCAGCCCCGTGTCAAGCAGGCCTTTGGCCCTGAAATCTGTATGTCCACCTTGAGGCC
 TGGGTAACCTTTGAGCAGCTGGTAACCTAGTTTACTCTTCATGATCCCAACCTCTTTGGTCCCCAG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001105790
Insert Size:	558 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001105790.1, NP_001099260.1</u>
RefSeq Size:	827 bp
RefSeq ORF:	558 bp
Locus ID:	287419
Cytogenetics:	10q24