

## Product datasheet for **RN204984**

### **Fosl2 (NM\_012954) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Fosl2 (NM\_012954) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Fosl2  
**Synonyms:** Fra2  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN204984 representing NM\_012954  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGTACCAGGATTATCCCGGAACTTTGACACCTCGTCCCGGGCAGCAGTGGCTCTCTGCGCACGCCG  
AGTCCTACTCCAGCGGTGGCGGGCCAGCAGAAGTTCGGGTAGATATGCCTGGCTCAGGCAGTGCCTT  
CATCCCCACAATCAACGCCATCACCACCAGCCAGGACCTGCAGTGGATGGTACAGCCCACAGTGATT  
ACCTCCATGTCCAACCCCTATCCGCGCTCGCACCCCTACAGCCCTGCCAGGCTTGCGATCAGTCCCTG  
GGCACATGGCTCTCCCCAGACCAGGAGTGATCAAGACCATTGGTACCACCGTGGGCCGAGAAGGAGAGA  
TGAGCAGCTGTCTCCGAAGAGGAGGAGAAGCGTCAATCCGGCGAGAGAGAAACAAGCTGGCCGCTGCC  
AAGTGTCCGAACCGTCGACGTGAGCTGACGGAGAAGCTGCAGACGGAGACAGAGGAGCTGGAGGAGAGA  
AGTCCGGTCTGCAGAAAGAAATCGCTGAGCTGCAGAAGGAGAAGGAGAAGCTGGAGTTCATGCTGGTGGC  
TCACGGCCCTGTGTGCAAAATCAGTCTGAGGAACGCCGATCTCCCCACCTCCGGGTACAGTCTTA  
CGCGGTACGGGCAGTGTGTCCGCCCTGTGGTGGTAAACAGGAGCCTCCTGAAGAGGACAGCCCCCTCT  
CCTCAGCAGGGATGGACAAGACCCAGCGCTGTGCATCAAGCCATTAGCATTGCTGGGGTGGTTCTTA  
CGGGGAAGAGCCTCTGCACACCCCATCGTGGTGACCTCCACACCTGCCATCACTCCGGGCACTTCAAAC  
CTTGCTTTCACCTACCCAGTGTCTGGAGCAGGAGTGCCTGCCTCGCCCTCCGAGTCTGTTCCAAGG  
CTCACCGCAGAAGCAGTAGCAGTGGGGACCAGTATCAGACTCCTTGAATCCCCACTCTGCTGGCTCT  
GTAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_012954



[View online »](#)

<b>Insert Size:</b>	984 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_012954.1</a> , <a href="#">NP_037086.1</a>
<b>RefSeq Size:</b>	984 bp
<b>RefSeq ORF:</b>	984 bp
<b>Locus ID:</b>	25446
<b>UniProt ID:</b>	<a href="#">P51145</a>
<b>Cytogenetics:</b>	6q14
<b>Gene Summary:</b>	member of the Fos family of transcription factors with circadian regulation [RGD, Feb 2006] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer protein (isoform 1).