

Product datasheet for **RN205015**

Scgb1d2 (NM_177482) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Scgb1d2 (NM_177482) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Scgb1d2
Synonyms:	Pbpc1bs; Psbpc1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN205015 representing NM_177482 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCACCATTAAGCTGAGCCTGTGTCTTCTGATCATGCTGGCTGTTTGTGCTATGAAGCTAATGCTA
GCCAGATCTGTGAACCTGTTGCCCATGAAACCATAAGCTTCTTAATGAAAAGTGAGGAAGAAGCTGAAGAA
GGAACCTGAGATGTATAATGCACCTCCAGCAGCTGTTGAAGCAAACTGGAAGTGAAGAGATGTGTAGAC
CAGATGAGCAATGGAGACAGATTGGTAGTAGCAGAAACTGGTATACATTTTTTTGGAATGTGGTGTGA
AACAAATGGGTAGAAACATATTATCCTGAGATCGATTTCTACTACGATATGAAC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_177482
Insert Size:	336 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_177482.1](#), [NP_803435.1](#)

RefSeq Size: 412 bp

RefSeq ORF: 336 bp

Locus ID: 309203

UniProt ID: [P02782](#)

Cytogenetics: 1q43

Gene Summary: member of the uteroglobin superfamily; major secretory glycoprotein of the rat ventral prostate gland [RGD, Feb 2006]