

Product datasheet for **RN205017**

Dtymk (NM_001106925) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dtymk (NM_001106925) Rat Untagged Clone
Tag: Tag Free
Symbol: Dtymk
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN205017 representing NM_001106925
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGC**C

ATGGCGTCGCGTCGCGGGGCGCTCATCGTGCTGGAGGGTGTGGACCGTGTGGCAAGACCACGCAGTGCC
TCAGGCTGGTGACCGCGCTGTGCGCCTCGGGCCACAGAGCGGAGCTGCTGCGTTCCCTGAAAGATCAAC
AGAAATCGGCAAGCTTCTGAGTTCCTACTTAGAAAAGAAAACCGAGCTCGAGGATCATTCCGTGCACCTG
CTCTTCTCTGCCAACCGCTGGGAACAAGTGCCATTAATTAAGGCGAAGTTGAACCAAGGGTGTGACCCCTG
TTTTGGACAGATACGCGTTTTCTGGGGTTCCTTCACTAGCGCCAAGGAGAATTTTTCTTGATTGGTG
CAAACAACCAGATGTGGGCCCTTCCCAAGCCTGATCTGATCCTGTTCTGACGTTACAATTGCTGGACGCT
GCTTCGAGGGGAGAGTTTGGTCTTGAGCGCTATGAGACTGGGACGTTCCAGGAGCAGGTTCTGCTGTGTT
TCCAGCAGCTCATGAGTGAGAAGAACCTCAACTGGAAGGTGGTTGATGCTTCCAAAAGCATTGAGGATGT
CCATGAAGAAATCCGTGCGCTCTCAGAGGACGCCATCCAAAATGCTGCACAGAGGCCACTGGAGGAGCTA
TGGAA**TAA**

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-RsrII
ACCN: NM_001106925
Insert Size: 639 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001106925.1](#), [NP_001100395.1](#)

RefSeq Size: 1235 bp

RefSeq ORF: 639 bp

Locus ID: 301622

Cytogenetics: 9q36