

## Product datasheet for RN205086

### Hspbp1 (NM\_139261) Rat Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Hspbp1 (NM_139261) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Hspbp1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN205086 representing NM_139261 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCAGACAAAGGCTCAGGGGGCAGTCGCCTCCCCCTCGCACTGCCTCCGGCTCCCAGGGTTGCTCGT  
CAGGGAGCAGTGGCTCCTCGCGGGGGGCTCAGGCAACCTCGCCTTCCACGGAACCTCCAAGGCCTGCT  
GCAGATGGCTATTACTGCGGGCTCTGAGGAACCAGACCCTCCTCCAGAACCCATGAGCGAGGAGAGACGC  
CAATGGCTGCAGGAAGCCATGTCAGTGCCTCCGAGGCCAGCGGAAGAGGTAGAGCAGATGAAGAAT  
GCCTCCGGGTCTGTCCAGGCCACCCCCAACTGCTGGTGAAGCTGAAGTGGCCACTGACCCAGCAGGA  
GGTGAAGGGGCACTAGAGCTGCTGGCAGACCTGTGCGAGAACATGGACAATGCGGCAGATTTCTGCCAG  
CTGTCAGGCATGCACCTGCTGGTGGTGCATACCTGGAGGCAGGAGCTGCAGGGCTGCGTTGGAGGGCAG  
CACAGCTCATCGGCACGTGCAGTCAGAATGTTGCAGCCATCCAGGAGCAGGTGTTGGGCTTGGTGCCCT  
GCGTAAGCTACTTCGGCTGCTCGACCGGACTCCTGCGACACGGTACGAGTCAATGCTCTCTTCGCCATC  
TCCTGTCTTGCCGAGAGCAGGAGGCTGGATTGCTGCAGTTCTCCGCCTGGATGGATTCTCAGTGTGA  
TGCGGGCCATGCAGCAGCAAGTGCAGAAGCTCAAGGTCAAGTCAAGTCAAGTCAAGTCAAGTCAAGTCAAGT  
GGCCACCCTGAGCACAAGGAACCTTTGCTCCATGGGGATGGTCCAGCAGCTGGTGGCCCTTGTGAGG  
ACAGAACACAGTCTTTCCATGAGCATGTGCTTGAGCCCTGTGCAGCCTTGTGACAGACTCCCTCAGG  
GTGTTTCGAGAATGCCGGGAGCCTGAGCTGGGCCTGGAGGAACTGCTCCGTCATCGCTGCCAGCTGCTGCA  
GCAGCATGAGGAGTACCAGGAGGAGCTGGAGTTCTGTGAAAAGCTGTTACAGACCTGTTTCTCCAGCCCT  
ACGGATGACAGCATGGATCGG**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_139261



<b>Insert Size:</b>	1074 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_139261.1</a> , <a href="#">NP_640354.1</a>
<b>RefSeq Size:</b>	1465 bp
<b>RefSeq ORF:</b>	1074 bp
<b>Locus ID:</b>	246146
<b>UniProt ID:</b>	<a href="#">Q6IMX7</a>
<b>Cytogenetics:</b>	1q12
<b>Gene Summary:</b>	inhibits chaperone activity by preventing ATP binding [RGD, Feb 2006]