

Product datasheet for **RN205091**

Ccdc130 (NM_001037644) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ccdc130 (NM_001037644) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Ccdc130
Synonyms:	RGD1307598
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN205091 representing NM_001037644 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGGGGAAAGGAAAGGCCAGAACAAAGTACTACCTCCGGAATTCACCCCGAAAGCATGGCTCTCTCA
ACCGATACCACACAGCCACCCACTTCGAGAACGGCTCGGAAGCTCTCACAAGGCATCCTTATTATCCG
ATTTGAAATGCCATACAACATCTGGTGTGACGGCTGTAAGAACCACATCGGCATGGGTGTGCGCTACAAT
GCAGAAAAGAAGAAGGTTGGCAATTACTACACGACTCCGATCTACAGGTTCCGGATGAAATGCCACCTCT
GTGTCAACTACATCGAGATGCAACAGACCCTGCCAACTGTGACTATGTCATCGTGAGTGGTGCCAGTCTG
CAAGGAGGAGCGTTGGGATATGGAAGACAATGAGCAGGTGCTGACCACAGAGCATGAGAAGAAAGAGAAA
CTGGAGACAGATGCCATGTTCCGCCTAGAGCATGGTGAGGCTGACCGCAGCACGCTAAAGAAGGCCCTTC
CCACCCTCAGCCACATCCAGGAGGCACAGAATGCCTGGAAGGATGACTTTGCTCTGAATAGCATGTACG
GAGGCACTTCGGGAGAAAAAGAAAGCCATGCAGGAGGAGGAAGAGAAGGACCAGGCCCTGCAGGCGAAG
GCCAACCTAGCCATCCCTCTTGTCAGAGTCCGAGGATGACCGCAGGCTGGCTGCCCTGCTAAGGTTGC
ACACCCTGGACTCTTATGAGGACAAACAGAGAATGAAGCGGACAGAGATCATCCACCGATCATGGTTCCT
CTCAGCCCAGGGACCCAGCACTAGCAGCAGCAAAGCTAGCACTGTCCTGAAGAAGCTGTGCCGGGGCCGC
AGACCACCCAGGCAGCGCAGGGGACCGGGGACCTGGGCATAGTGAGGAGAAAGTCTCGAGAAGCTC
CAGAAAGTCCCCAGTGTACAGCGGACAACCTCTTGTGAGAAGAACCAAGAGGACCACCCGGAACCCCCC
AGACAGCAAGACACTTCAAGGGACAGCAGAGGCACCCAGAACAGCAAGACACTAGAATCTAAGAGAAAC
TGTAAGTATCAAGCCTTACCCCTAGGTTCTCCAGGAGGACCTGCTGCACCCCAACACCCCAACGCCT
CCCTGGTGGCTGACTACTCTGACTCAGAGAGTGAG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001037644
Insert Size:	1158 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001037644.1, NP_001032733.1</u>
RefSeq Size:	1590 bp
RefSeq ORF:	1158 bp
Locus ID:	304656
UniProt ID:	<u>Q32PZ9</u>
Cytogenetics:	19q11