

Product datasheet for **RN205100**

Vom1r7 (NM_001008907) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Vom1r7 (NM_001008907) Rat Untagged Clone

Tag: Tag Free

Symbol: Vom1r7

Synonyms: V1re3

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN205100 representing NM_001008907
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGATTTTAGGGATATGGTAATAGGAATACTGCTTCACTTGAGATTTTACTTGGTATTATTGGAACT
TCTCCCTCTTTCTACTATCTAATCCTTTACTATAAGGAACACACATTAATAATCATAGATACAATTCT
TATACATGTGTTTCATATCTAATCCTTGATCCTTCTCTCAAAGGACCCCTGAGGTACTGGGAGCTTTT
GGGTCAAAGCAGTTATTCAATGATGTTGAGTGCAAAATTATCTTATATGTTCAAAGAATTGCCAGATGTA
TGTCATAAGCACCACTTGCTCTTGAGTGTCTTCAAGCAATCACCATCAGCCCGAGGAAGTCTTGCTG
GAAGGAATTTAAAGTCAAACCACAAAGTTATGGGACTGTCTATTTCCCTCTGCTGGATTATGTACATG
TTACTAAATATGCTTTTCCCTTTATATACATCTACCAAAGGTAATAGAGAAAATCGGACCCAAAAGAGTG
ATTCTGAATTCTGCTACTCTGGAGGCCGTGACAAAATAGTAGATTTATTGTATGCAGTATTTTGTGTGT
TCCAGAAGTCTTATTTTCTACTCATTGTAGTTCCAGTACCTTCATGATTGTTATTCTTTATAGACAC
AAGAAAAGGGTTCAACACATCCTTCACATTCATGCTTCTTTAGAATCTCTGCTGAAAACAGAGCCACAC
AGACCATTTTGGTCTTGGTTTTACCTTTCTAGTGTTTTATACTCTTTCTCAATTTTACAAGGCTATAT
TGCTCTTTACATGATCCCAATTGGTGGTAATGAATATCACAGCTATCATAGCTTTGTGTTTCCCTACT
TTAGGCCCATTTGTGATGAGTCATAATTTTACTGTTTCCGGATTTTGCTTTACCTGGATAAAAAACATTT
AA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001008907



[View online »](#)

Insert Size:	912 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001008907.1</u> , <u>NP_001008907.1</u>
RefSeq Size:	912 bp
RefSeq ORF:	912 bp
Locus ID:	494259
Cytogenetics:	1q12