

Product datasheet for RN205152

Abhd4 (NM_001108866) Rat Untagged Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Abhd4 (NM_001108866) Rat Untagged Clone |
| Tag: | Tag Free |
| Symbol: | Abhd4 |
| Synonyms: | Abh4 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >RN205152 representing NM_001108866 Red=Cloning site Blue=ORF Orange=Stop codon |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGCTGGCTCAGCTCGACCCGGCAGGGCTTGTACTATGGCTGATGATCTGGAGCAGCAGCCCCAGG
GCTGGCTGAGTAGCTGGCTCCCCACTTGGCGCCCCACTTCCATGTCTCAGCTGAAGAATGTGGAAGCCAG
GATCCTCCAGTGTCTCCAGAACAAGTTCTGGCCGTTATGTATCCCTCCAAACCAGAACAAGATCTGG
ACGGTACTGTGAGCCCAGAGCAAAAGGACCGCACCCCTCTGGTATGGTACATGGCTTTGGGGTGGCG
TGGCCTCTGGATCCTCAACATGGATCCCTGAGTGCCCGCCGACACTCCATACCTTTGATCTGCTTGG
TTTTGGGCGAAGCTCAAGGCCTACGTTCCCAAGGACCCAGAAGGAGCTGAAGATGAGTTTGTGACCTCA
ATAGAGACATGGCGGGAGACAATGGGAATCCCTACCATGATCCTCTGGGCGACAGTTTGGGAGGATTCC
TGGCCACTTCTACTCTAACAAGTACCCTGAAAGAGTTAAACATCTAATCCTGGTGGATCCCTGGGGCTT
TCCCTTACGACCGACTGACCCAGTGAGATCCGTACACCTCCAACCTGGGTCAAGGCTGTGGCATCTGTC
CTGGGACGCTCCAATCCACTGGCTGTTCTCGAGTGGCTGGGCCCTGGGGCCTGGGCTGGTGCAGAGAT
TCCGTCCAGACTCAAGCGCAAATTTGCAGACTCTTTGAAGATGATACCATCTCGGAATACATCTACCA
CTGCAATGCACAGAATCCAGTGGGAAACGGCATTCAAAGCCATGATGGAGTCCCTCGGCTGGGCCCGG
CGCCCCATGTTGGAGCGAATCCACTTGATTTCGAAAGATGTGCCCATCACCATGATCTATGGGGCCAATA
CCTGGATCGACACCAGCACAGGGAAAAAGGTGAAGCTGCAAAAGGCCAGACTCCTACGTCGAGACATGGA
GATCGAGGGTGCATCGACCATGTCTACGCTGACCAGCCACACATCTTCAATGCTGTGGTGGAGGAGATC
TGTAACTCAGTTGACTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online >](#)

| | |
|-------------------------------|---|
| ACCN: | NM_001108866 |
| Insert Size: | 1068 bp |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001108866.1</u> , <u>NP_001102336.1</u> |
| RefSeq Size: | 2399 bp |
| RefSeq ORF: | 1068 bp |
| Locus ID: | 364380 |
| Cytogenetics: | 15p13 |