

## Product datasheet for RN205178

### Kdelr1 (NM\_001017385) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Kdelr1 (NM\_001017385) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Kdelr1  
**Synonyms:** MGC109169  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >RN205178 representing NM\_001017385  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAACCTCTCCGATTCCTGGGAGACCTCTCCACCTCCTAGCTATCATCTTGTACTGCTCAAATCT  
GGAAGTCCCGATCGTGTGCGGGCATTTCAGGGAAGAGCCAAGTCTGTTGCTGTGGTGTCACTGCCCG  
ATATCTGGACCTGTTACCAACTACATTTCACTCTACAATACCTGCATGAAGGTGGTCTACATAGCTTGT  
TCCTTTACCACGGTGTGGATGATCTACAGCAAGTCAAAGCTACGTATGATGGGAACCACGACACTTTCC  
GGGTGGAGTTCCTCGTTTCTACTGCTGTTTTGGCGTTCCTGGTTAACCATGACTTCACCCCTCTGGA  
GATCCTCTGGACTTTCTCCATCTACCTGGAGTCACTGGCTATCTTGCCACAGCTGTTATGGTGAGCAAG  
ACGGGCGAGGCGAGACCATCACCAGCCATTACCTGTTCCGACTGGGTGTCTACCGTACACTCTATCTCT  
TCAACTGGATCTGGCGCTACCACTTCGAGGGCTTCTTTGACCTCATCGCCATTGTTGCTGGCCTGGTCCA  
GACAGTCTCTACTGCGATTTCTTCTACCTCTACATCACCAAAGTCTAAAGGGGAAGAAGCTGAGCTTG  
CCGGCG**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001017385  
**Insert Size:** 639 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001017385.1</a></u> , <u><a href="#">NP_001017385.1</a></u>
<b>RefSeq Size:</b>	1476 bp
<b>RefSeq ORF:</b>	639 bp
<b>Locus ID:</b>	361577
<b>UniProt ID:</b>	<u><a href="#">Q569A6</a></u>
<b>Cytogenetics:</b>	1q22
<b>Gene Summary:</b>	Receptor for the C-terminal sequence motif K-D-E-L that is present on endoplasmic reticulum resident proteins and that mediates their recycling from the Golgi back to the endoplasmic reticulum.[UniProtKB/Swiss-Prot Function]