

## Product datasheet for **RN205181**

### Mrgprd (NM\_001001506) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Mrgprd (NM\_001001506) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Mrgprd  
**Synonyms:** Mgrd; Mrgd; TGR7  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >RN205181 representing NM\_001001506  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAACACTCCTTATAGCAGCCCAGCCCCAGGTCTGACCATCAGCCCCACCATGGACCCTGTGACCT  
GGGTTTACTTTTCAGTGACATTCCTGGCCATGGCCACCTGTGTGTGGGATAGTGGCAACTCCATGGT  
GATTTGGCTACTGAGTTCCACAGGGTGCAGAGGTCCCCCTTCGACCTACGTGCTCAACTGGCGGTG  
GCCGACCTCCTCTCCTGCTCTGCATGGCCTCCCTGCTCAGTCTGGAAACAGGGCCCTGCTCACAGCCA  
GCACCTCCGCCAGAGTCTACGAGGGGATGAAGAGAATCAAGTACTTTGCCTACACAGCAGGCCTGAGCCT  
GCTGACGGCCATCAGCACCCAGCGCTGTCTCCTGGTCTTTCCCATCTGGTATAAGTGCCACCGGCC  
CAGCACCTGTCGGGGTGGTATGTGGTGTGCTGTGGGCACTGGCCCTCCTGATGAACTTCTGGCTTCT  
TCTTCTGTGTTCAATTCTGGCATCCCGACAAATACCAGTGCTTCAAGGTGGACATGGTTTTCAACAGTCT  
TATCCTGGGGATCTTCATGCCCGTCATGGTCTGACCAGCGCCATCATCTTCATCCGCATGCGAAAGAAC  
AGCCTGCTGCAGAGACGGCAGCCTCGGCGGCTCTACGTGGTCACTGACTTCCGTCTTGTCTTCTTA  
CCTGTTCTCTGCCGTTGGGCATCAACTGGTCTTACTCTACTGGTGGAAGTCCCGCAGGCCGCTGAGGCT  
CCTGTACGTCTGCTCATCAGCTTCTCCTCGTCTTTGAGCAGCAGGCCAACCCAGTCACTACTTCCCTC  
GTGGGCAGCCAGAAGAGCCACCGGCTGCAGGAGTCTCTGGGTGCTGTGCTGGGGCGGCACCTTCAGGACG  
AGCCTGAAGGCAGGGAGACGCCATCCACATGTAATAATGATGGGGT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Chromatograms:** [https://cdn.origene.com/chromatograms/ja2520\\_e01.zip](https://cdn.origene.com/chromatograms/ja2520_e01.zip)

**Restriction Sites:** SgfI-MluI



[View online >](#)

<b>ACCN:</b>	NM_001001506
<b>Insert Size:</b>	960 bp
<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_001001506.1</a> , <a href="#">NP_001001506.1</a>
<b>RefSeq Size:</b>	960 bp
<b>RefSeq ORF:</b>	960 bp
<b>Locus ID:</b>	293648
<b>UniProt ID:</b>	<a href="#">Q7TN41</a>
<b>Cytogenetics:</b>	1q42
<b>Gene Summary:</b>	seven transmembrane domain G-protein coupled receptor; most often found in sensory neurons [RGD, Feb 2006]