

Product datasheet for **RN205247**

Scml4 (NM_001107638) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Scml4 (NM_001107638) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Scml4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN205247 representing NM_001107638 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCCCTCACACACACTTTCAAAGTGTCTGGACCTGTTATTGAAATAAGCCTGGCACCAGCAGAAATGC
 AGACCCAGAGGGTCCCAGGGAGGAAGCGTGGCCGGCTCCACTTCACTCCACGCGGGTGCAGATGGCAGT
 GCATAACCTTTATTCTGCTTCTTCTGCCTCCGTGCCAGCAGTGACAATCCCGAAGAAGAGAGGACGGAAA
 CCCAGGTACAAGATCAAGTCTCCGGTGTCTGACTCCCTTGGCCCTGTCACTCCACGGAGCACTCCAG
 AACCCGACCTCAGCTCCATCCCTCAGGATGCAGCACTGTCCCGAGCTTGGTGGTCCCACAGGCTGTCAC
 AGTCTCTGCTTCCCTTTGATGGAAGCAGCACCTGAGGAATCTGCCTGTGGTGAACAGTGTGGGCTACGTC
 CTCCGCTTCTGACCAAACCTGTGCCGAGCCTCCTGTGTGACAACCTCTTCAGCCACCTGCCTTTTCCTG
 GGAGCATCGGTGCCTCTGAGAAGGCCAGGAGAAAGAGGATGGGAGGACTGAGTCAGCCAAGGTAGCCAC
 TGCTGAAGAGTGTCTGGTAAACGCAGTGGGAATGAACCGCTACACCATAGATGTCTCTCCACTGCCTTT
 AGCCACAGGGGCTCTGTGACCACTCTTCTCGCTGTACTACAAGAGGCTGAACCTCTGGAGACAGTCACC
 TCGTGGGAGGCCCGCCACTACCACAGTGGTTCCCGTACCAACCCTGTGCCCTCTGGAGGTTTCATCAAC
 ACCTGGGCTGAGGCTCCAGCTTCCAGCCCCAAGAGGAATGGTACAGCCATCGAAGGAAACAGATGTGCC
 CCGAGCCCATCTCCAGAGATACAGGATACCAGGAGGCCAAGCAGCAGGAATCCTTCTACCTGGACCGTGG
 AGGACGTGGTCCGGTTTGTGAAAGATGCAGACCTCAGGCCCTGGGGCTCATGTGGAGCTTTCAGAAA
 GCATGAGATTGATGGCAATGCCTTGTGTTGCTAAGGAGTGACATGATCATGAAATACCTGGGACTGAAG
 CTGGGCCCCGCGCTGAACTCTGCTATCACATTGACAAGCTGAAGCAAGCTAAGTTC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001107638


[View online »](#)

Insert Size:	1110 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001107638.1, NP_001101108.1</u>
RefSeq Size:	4333 bp
RefSeq ORF:	1110 bp
Locus ID:	309859
Cytogenetics:	20q13