

Product datasheet for **RN205256**

Psmid7 (NM_001107426) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Psmid7 (NM_001107426) Rat Untagged Clone
Tag: Tag Free
Symbol: Psmid7
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN205256 representing NM_001107426
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCCGGAGCTGGCGGTGCAGAAAGTGGTGGTCCACCCCTGGTGTCTCAGTGTGGTGGATCATTCA
ACCGAATCGGCAAGTTGGAAACCAGAAGCGGGTTGTTGGTGTGCTTTTGGGATCATGGCAAAGAAAGT
CCTTGATGTATCCAACAGCTTTGCAGTACCTTTTGATGAAGATGACAAAGATGATTCTGTCTGGTTTTTA
GACCATGATTATTTGAAAACATGTATGGGATGTTTAAGAAGGTCAATGCCAGAGAAAGAATAGTTGGGT
GGTACCACACAGGCCCAAACACTGCAAGAACGATATCGCCATCAATGAAGTCAAGAGATACTGCC
CAACTCGTACTGGTCAATATTGATGTGAAGCCAAAAGACTTGGGGCTCCGACCGAAGCCTACATCTCA
GTGGAGGAAGTTCATGATGACGGAACACCAACATCAAAAACCTTTGAGCACGTGACTAGTGAAATTGGAG
CAGAGGAGGTGAGGAGTTGGAGTGGAGCACTTGTACGGGACATCAAGGACACTACAGTGGGGACTCT
GTCTCAGAGGATCACTAACCAGTCCATGGTTTGAAGGGACTGAACTCCAAGCTCCTGGATATCAGGACC
TACCTGGAGAAGGTAGCCAGCGCAAGCTACCCATCAACCACCACATCATATACCAGCTGCAGGACGTCT
TCAACCTGCTGCCGGACGCCAGCCTGCAGGAGTTCGTCAAGGCCCTTCTACCTGAAGACCAATGACCAGAT
GGTTGTGGTGTACCTGGCCTCGCTGATCCGCTCTGTGGTGCCTTGATACCTCATCAACAACAAGATT
GCCAACCGCATGCCGAGAGAAAGAGGGACAGGAGAAGGAGGAGAGCAAGAAGGAGAGAAAAGACGACA
AGGAGAAGGAGAAGAGCGATGCCAAGAAAGAGGAGAAAAAGGAGAAAAAG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001107426
Insert Size: 963 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001107426.1</u> , <u>NP_001100896.1</u>
RefSeq Size:	1588 bp
RefSeq ORF:	963 bp
Locus ID:	307821
Cytogenetics:	19q12