

Product datasheet for **RN205276**

Tgfbr2 (NM_031132) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tgfbr2 (NM_031132) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Tgfbr2
Synonyms:	TGF-beta 2; Tgfbr2T
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >RN205276 representing NM_031132
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGTCGGGGCTGCTCCGGGCTGTGGCCGCTGCACATCGTCTGTGGACGCGCATCGCCAGCACGA
 TCCCGCCGCACGTTCCCAAGTCGGTTAACAGCGATCTGATGGCGGGTGACAAACAGCGGTGCGGTCAAGCT
 TCCACAGCTGTGCAAGTTTTGCGACGTGACACTGTCCACTTGTGACAACCAGAAGTCTTGATGAGCAAC
 TGCAGCGTCACCTCCATCTGTGAGAAGCCGAGGAAGTCTGCGTGGCCGTGTGGAGGAAGAACGACAAGA
 ACATTACTCTGGAGACGTTTGGCACGACCCCAAGTTCACCTACCACGGCTTCACTCTGGAAGATGCCAC
 TTCTCCCACGTGCGTCATGAAGGAAAAGAAAAGGGCGGGCGAGACCTTCTTCATGTGCTCTGTAACACA
 GAGGAGTGAACGATTACATCATCTTTAATGAAGAATACACCACCAGCAGTCTGACCTGTTGCTGGTCA
 TTATCAAAGTGACGGGCGTCAGCCTCCTGCCTCCGCTGGGGATTGCCATAGCTGTCATTGCCATCTTCTA
 CTGTTACCGTGTCCACCGCAGCAGAAGTTGAGCCCTCCTGGGAGAGCAGCAAGCCCCGGAAGCTCATG
 GATTTACGCGACAACGCGCCATCATCCTGGAGGACGACCGCTCTGACATCAGCTCTACGTGCGCCAACA
 ACATCAACCACAATACGGAAGTCTGCCATCGAGCTGGACACGCTGGTGGGGAAGGGCCGCTTCGCCGA
 GGTCTACAAGGCCAAGCTGAAGCAGAACACTTCAGAGCAGTTTGGAGCCGTGGCTGTCAAGATCTTCCCC
 TACGAGGAATACTCCTCGTGGAAAACGGAGAAGGACATCTTCTCTGACATCAACCTGAAACACGAGAACA
 TCCTGCAGTTCCTGACGGCCGAGGAGCGGAAGACGGAGATGGGCAAGCAGTACTGGCTGATCACGGCGTT
 TCATGCTAAGGGCAACCTGCAGGAGTACCTCACTAGGCACGTATCAGCTGGGAGGACCTGCGGAAGCTG
 GGCAGCTCCCTGGCCCGGGCATCGCCATCTGCACAGTGACCACACTCCTTGTGGGAGGCCAAGATGC
 CCATTGTTACAGGGACCTCAAGAGCTTAACATCCTAGTAAGAAGCGATCTAACCTGTTGCCTGTGTGA
 CTTCCGGTTGTCCTTGGCCCTGAGCCCTACTCTGTCTGTGGATGACCTGGCCAACAGCGGGCAGGTGGGA
 ACAGCGAGATACATGGCCCCGGAAGTCTAGAATCCAGGATGAATTTGGAGAACATGGAGTCTTCAAGC
 AGACGGATGTCTACTCCATGGCTCTAGTACTCTGGGAAATGACGTGCGCTGCAACGCTGTGGGAGAAGT
 GAAAGATTATGAGCCCCGTTTGGTTCCAGAGTGCGGGAGCACCCCTGTGTGGAGAGCATGAAAGACAAC
 GTGCTGAGAGACCGAGGACGGCCTGAAATCCCCAGCTTCTGGCTCAACCACCAGGGCATCCAGATCGTGT
 GTGAGACTGACCGAGTGTGGACCACGACCCCGAGGCCGCTCTCACAGCCAGTGTGTGGCAGAGCG
 CTTCACTGAACTGGAGCACCCAGATAGACTCTCCGGGAGGAGCTGCTCCAGGAGAAGATTCCAGAAGAC
 GGCTCCCTGAACACTACCAAA**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAAAACATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja3128_a12.zip

Restriction Sites: SgfI-MluI

ACCN: NM_031132

Insert Size: 1704 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_031132.3](#), [NP_112394.3](#)

RefSeq Size: 2080 bp

RefSeq ORF: 1704 bp

Locus ID: 81810

UniProt ID: [P38438](#)

Cytogenetics: 8q32

Gene Summary: receptor for TGF-beta; may play a role in lung maturation and response to hypoxia [RGD, Feb 2006]